

# Archived Document

This archived document is no longer being reviewed through the CLSI Consensus Document Development Process. However, this document is technically valid as of January 2017. Because of its value to the laboratory community, it is being retained in CLSI's library.



April 2008

## LIS01-A2

### Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition

This document describes the electronic transmission of digital information between clinical laboratory instruments and computer systems.

A standard for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement, but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advancements in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential, and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeals Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeals, documented in the CLSI Standards Development Policies and Processes, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute  
950 West Valley Road, Suite 2500  
Wayne, PA 19087 USA  
P: +1.610.688.0100  
F: +1.610.688.0700  
[www.clsi.org](http://www.clsi.org)  
[standard@clsi.org](mailto:standard@clsi.org)

ISBN 1-56238-665-4  
ISSN 0273-3099

LIS01-A2  
Vol. 28 No. 13  
Replaces LIS1-A  
Vol. 23 No. 7

---

# Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition

Volume 28 Number 13

David Chou, MD  
Andrzej J. Knafel, PhD  
Charles D. Hawker, PhD, MBA, FACB  
Ed Heierman, PhD  
David A. Lacher, MD, MEd  
William Neeley, MD, FACP  
Eugene T. Reilly  
Richard S. Seaberg, MT(ASCP)

## Abstract

CLSI document LIS01-A2—*Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition* describes the electronic transmission of digital information between clinical laboratory instruments (those that measure one or more parameters from one or multiple samples) and computer systems (those that are configured to accept instrument results for further processing, storage, reporting, or manipulation).

Clinical and Laboratory Standards Institute (CLSI). *Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition*. CLSI document LIS01-A2 (ISBN 1-56238-665-4). Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2008.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org). If your organization is not a member and would like to become one, and to request a copy of the catalog, contact us at: Telephone: 610.688.0100; Fax: 610.688.0700; E-Mail: [customerservice@cls.org](mailto:customerservice@cls.org); Website: [www.clsi.org](http://www.clsi.org).



CLINICAL AND  
LABORATORY  
STANDARDS  
INSTITUTE®

Copyright ©2008 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, companion product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedure manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

### **Suggested Citation**

CLSI. *Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition*. CLSI document LIS01-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.

### **Previous Edition:**

April 2003

### **Reaffirmed:**

September 2011

### **Archived:**

January 2017

ISBN 1-56238-665-4

ISSN 0273-3099

## Committee Membership

### Area Committee on Automation and Informatics

**David Chou, MD**  
**Chairholder**  
 University of Washington Medical  
 Center  
 Seattle, Washington

**Andrzej J. Knafel, PhD**  
**Vice-Chairholder**  
 Roche Instrument Center AG  
 Rotkreuz, Switzerland

Charles D. Hawker, PhD, MBA,  
 FACB  
 ARUP Laboratories, Inc.  
 Salt Lake City, Utah

Ed Heierman, PhD  
 Abbott Diagnostics  
 Irving, Texas

David A. Lacher, MD, MEd  
 National Center for Health Statistics  
 Hyattsville, Maryland

William Neeley, MD, FACP,  
 DABCC  
 Detroit Medical Center University  
 Laboratories  
 Detroit, MI

Eugene T. Reilly  
 FDA Ctr. For Devices/Rad Health  
 Rockville, MD

Richard S. Seaberg, MT(ASCP)  
 North Shore University Hospital  
 Manhasset, New York

#### Advisors

Michael G. Bissell, MD, PhD, MPH  
 Ohio State University  
 Columbus, Ohio

Suzanne H. Butch, MA,  
 MT(ASCP), SB  
 The University of Michigan  
 Ann Arbor, Michigan

Randy R. Davis  
 Siemens Medical Solutions  
 Diagnostics  
 Bear, Delaware

Al DeStefano  
 Sysmex America, Inc.  
 Mundelein, Illinois

Jeffrey A. DuBois, PhD  
 NOVA Biomedical Corp.  
 Waltham, Massachusetts

Arden W. Forrey, Jr., PhD, FACB  
 University of Washington  
 Seattle, Washington

David A. Herold, MD, PhD  
 VA (San Diego) Medical Center  
 San Diego, California

Georg E. Hoffmann, MD  
 Trillium GmbH  
 Grafrath, Germany

Stephen Howlett  
 Beckman Coulter, Inc.  
 Miami, Florida

Brian Richard Jackson, MD  
 ARUP Laboratories  
 Salt Lake City, Utah

Gary W. Kramer, PhD  
 National Institute of Standards and  
 Technology  
 Gaithersburg, Maryland

Rodney S. Markin, MD, PhD  
 Univ. of Nebraska Medical Center  
 Omaha, Nebraska

Richard A. McPherson, MD  
 Medical College of Virginia  
 Hospital  
 Richmond, Virginia

Paul J. Mountain, MSc, MT(ASCP)  
 Flamborough, Ontario, Canada

David O'Bryan, PhD  
 Hibernia Consulting  
 Portsmouth, New Hampshire

Jeff Quint, PhD  
 Beckman Coulter, Inc.  
 Brea, California

Hiroshi Sekiya  
 Olympus America Inc.  
 Irving, Texas

Russell H. Tomar, MD  
 Chicago, Illinois

#### Staff

Clinical and Laboratory Standards  
 Institute  
 Wayne, Pennsylvania

Lois M. Schmidt, DA  
*Vice President, Standards  
 Development and Marketing*

David E. Sterry, MT(ASCP)  
*Staff Liaison*

Patrice E. Polgar  
*Project Manager*

Melissa A. Lewis  
*Editor*

### Acknowledgment

CLSI and the Area Committee on Automation and Informatics gratefully acknowledge the contributions of Ed Heierman, PhD, Abbott Diagnostics; Randy Davis, Siemens Medical Solutions Diagnostics; and Bill Coughlin, BS, MBA, Data Innovations, for their contributions during the revision of the LIS01-A2 consensus standard.



## Contents

Abstract.....	i
Committee Membership.....	iii
Foreword.....	vii
1 Scope.....	1
2 Introduction.....	1
3 Terminology.....	1
4 Significance and Use.....	2
5 Physical Layer for Serial Binary Data Exchange.....	3
5.1 Overview.....	3
5.2 Electrical Characteristics .....	3
5.3 Mechanical Characteristics .....	5
6 Data Link Layer for Serial Binary Data Exchange .....	5
6.1 Overview.....	5
6.2 Establishment Phase (Link Connection).....	6
6.3 Transfer Phase.....	7
6.4 Termination Phase (Link Release).....	9
6.5 Error Recovery.....	9
6.6 Restricted Message Characters .....	10
7 Physical Layer for TCP/IP Data Exchange.....	10
7.1 Overview.....	10
7.2 Electrical Characteristics .....	10
7.3 Mechanical Characteristics .....	12
8 Data Link Layer for TCP/IP Data Exchange .....	12
8.1 Overview.....	12
8.2 Establishment Phase (Link Connection).....	13
8.3 Transfer Phase.....	14
8.4 Termination Phase (Link Release).....	16
8.5 Error Recovery.....	16
8.6 Restricted Message Characters .....	17
References.....	18
Appendix A. Mandatory Information .....	19
Appendix B. Nonmandatory Information .....	20
Summary of Delegate Comments and Area Committee Responses .....	21
The Quality Management System Approach.....	24
Related CLSI Reference Materials .....	25





## Foreword

In 2001, ASTM Committee E31 decided to restructure its operations, with the intent of focusing on standards-development issues such as security, privacy, and the electronic health record. Part of the reorganization plan was to transfer responsibility for E31.13 standards to CLSI, then known as NCCLS.

Following this transfer, nine standards (formerly ASTM E792; E1029; E1238; E1246; E1381; E1394; E1466; E1639; and E2118) were redesignated as CLSI/NCCLS standards LIS1 through LIS9.<sup>1-8</sup> This collection of standards provides a wide variety of information relating to clinical laboratory computer systems. Some included documents are of general interest as reference sources; others represent specifications of primary importance to instrument manufacturers. LIS2 is a revision of the former ASTM E1381-02.

The Area Committee on Automation and Informatics has assumed responsibility for maintaining the documents and will revise or update each document in accord with the CLSI Administrative Procedures. The area committee prioritized LIS1-A as the second standard from this collection to be updated, incorporated into the CLSI document template, and advanced through the CLSI consensus process. The area committee will revise other documents in the series in a similar manner.

With the transfer of the former ASTM standards, the Area Committee on Automation and Informatics has expanded its mission statement to include laboratory information systems. In the future, the area committee will develop additional standards addressing informatics issues, as well as issues related to the integration of patient clinical data.

This document replaces the first edition of the approved guideline, LIS1-A, which was published in 2003. Several changes were made in this edition; among them, TCP/IP communication is now included (Sections 4.4 and 4.5) and the state diagram was replaced (see Appendix A) so it is consistent with the text of the document.

The revisions in this edition of the LIS01 standard are also intended to delineate this document from its former ASTM edition. The title and text have been revised throughout to indicate that this standard applies to clinical laboratory instruments. The term *computer* has been replaced with the term *information* to better reflect the current terminology (ie, LIS) and the headings of Sections 6 and 8 have been changed to make them more specific.

## Key Words

data link layer, physical layer, serial binary data exchange, TCP/IP data exchange



# Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard— Second Edition

## 1 Scope

This specification describes the electronic transmission of digital information between clinical laboratory instruments and computer systems.

This specification addresses the low-level protocol used for both serial binary data exchange and TCP/IP data exchange. For message content in the interface between clinical instruments and computer systems, reference CLSI/NCCLS document LIS2.<sup>1</sup>

## 2 Introduction

The clinical laboratory instruments under consideration are those that measure one or more parameters from one or more patient samples. Often they will be automated instruments that measure many parameters from many patient samples. The computer systems considered here are those that are configured to accept instrument results for further processing, storage, reporting, or manipulation. This instrument output may include patient results, quality control results, and other related information. Typically, the computer system will be a clinical laboratory information management system (CLIMS).

The terminology of the International Organization for Standardization (ISO) Reference Model for Open Systems Interconnection (OSI) is generally followed in describing the communications protocol and services.<sup>9</sup> The electrical and mechanical connection between instrument and computer is described in the Physical Layer sections (see Sections 5 and 7). The methods for establishing communication, error detection, error recovery, and sending and receiving of messages are described in the Data Link Layer sections (see Sections 6 and 8). The data link layer interacts with higher layers in terms of sending and receiving “messages,” handles data link connection and release requests, and reports the data link status.

## 3 Terminology

**3.1 receiver** – the device that responds to the sender and accepts the message.

**3.2 sender** – the device that has a message to send and initiates the transmission process.

**3.3** The parts of a communication between instrument and computer are identified by the following terms. The parts are hierarchical and are listed in order of most encompassing first.

**3.3.1 session** – a total unit of communication activity, used in this standard to indicate the events starting with the establishment phase and ending with the termination phase, as described in subsequent sections.

**3.3.2 message** – a collection of related information on a single topic, used here to mean all the identity, tests, and comments sent at one time; **NOTE:** When used with CLSI/NCCLS document LIS2,<sup>1</sup> this term means a record as defined by CLSI/NCCLS document LIS2.<sup>1</sup>

**3.3.3 frame** – a subdivision of a message, used to allow periodic communication housekeeping, such as error checks and acknowledgments.

## 4 Significance and Use

- 4.1** Nearly all recent major clinical instruments have provisions for connection to a computer system, and in nearly all laboratories that have implemented a CLIMS, there is a need to connect the laboratory's high volume automated instruments to the CLIMS, so results can be transferred automatically. To accomplish this connection, both the instrument and the computer must have compatible circuits and appropriate software, and there must be a proper cable to connect the two systems.
- 4.1.1** Without this standard specification, the interface between each different instrument and each different computer system is likely to be a different product. This increases the cost, the chances for compatibility problems, and the difficulty of specifying and designing a proper system. In addition, interfaces for every instrument-computer combination may not be available, forcing expensive and time-consuming custom development projects.
- 4.2** This standard specification defines the electrical parameters, cabling, data codes, transmission protocol, and error recovery for the information that passes between the instrument and the laboratory computer. It is expected that future products from instrument manufacturers and computer system developers, released after the publication of this specification, will conform to this specification, and that will lead to plug-together compatibility of clinical instruments and computer systems.
- 4.3** This standard specification does not define the transmission of character encodings that require greater than eight bits; however, character encodings greater than eight bits can be transmitted through the use of appropriate algorithmic transformations. Examples include using Shift JIS for the transmission of Japanese characters defined by the Japanese Industrial Standards JIS X 0207 (single-byte characters) and JIS X 0208 (double-byte characters), as well as the transmission of the Unicode/UCS-2 character set with a UTF-8 encoding. These techniques do not violate the character-oriented protocol and associated character restrictions described by the data link layer of this standard. There are no provisions in the standard to identify if a specific transmission approach is being used; therefore, the character encoding and transmission approach must be agreed upon outside of the context of this standard.
- 4.4** This standard specification was originally intended for serial communication between clinical laboratory instruments and computer systems. In order to address bandwidth concerns with serial communication and to address a shift toward the use of networked solution, this standard was updated to include TCP/IP communication as well as frame sizes of 64 000 characters.
- 4.4.1** It is recommended that the use of large frame sizes be restricted to TCP/IP communication, where the frames will be broken down into smaller packets for transmission. This will ensure proper error detection and recovery for the transmissions.
- 4.4.2** When the TCP/IP data exchange is used, additional overhead is encountered, because TCP/IP is itself a low-level protocol; and LIS01-A2 protocol activities such as error checking and error recovery are also performed as part of the TCP/IP protocol. In addition, this standard specification is often used in combination with CLSI/NCCLS document LIS2,<sup>1</sup> which defines specific message formats. It may be more appropriate to use CLSI/NCCLS document LIS2<sup>1</sup> directly over TCP/IP, the Health Level Seven Protocol for electronic exchange, or the CLSI/NCCLS standard AUTO3<sup>10</sup> for communication with automated clinical laboratory systems, instruments, devices, and information systems when using the TCP/IP communication protocol to exchange information.
- 4.5** Although this specification addresses serial communication using a 25-pin connector and TCP/IP communication, it can be adapted to other standard communication approaches such as a 9-pin

serial connector, wireless, high speed USB, and firewire. The physical layer specifications for other approaches must be agreed upon outside the context of this standard. Adaptation to 9-pin, USB, etc. is standard, and can be accomplished in a manner transparent to the specification.

## 5 Physical Layer for Serial Binary Data Exchange

### 5.1 Overview

The mechanical and electrical connection for serial binary data bit transmission between instrument and computer system is described in the physical layer. The topology is point-to-point, a direct connection between two devices.

### 5.2 Electrical Characteristics

The voltage and impedance levels for the generator and receiver circuits are as specified in the EIA-232-D-1986 standard.<sup>11</sup>

#### 5.2.1 Signal Levels

5.2.1.1 For the data interchange circuits, a marking condition corresponds to a voltage more negative than minus three volts with respect to signal ground at the interface point. A spacing condition corresponds to a voltage more positive than plus three volts with respect to signal ground at the interface point.

5.2.1.2 Binary state ONE (1) corresponds to the marking condition; binary state ZERO (0) corresponds to the spacing condition.

5.2.1.3 The signal levels conform to the EIA-232-D-1986 standard.<sup>11</sup>

#### 5.2.2 Character Structure

5.2.2.1 The method of data transmission is serial-by-bit start/stop. The order of the bits in a character is:

- (1) One start bit, corresponding to a binary 0;
- (2) The data bits of the character, least significant bit transmitted first;
- (3) Parity bit; and
- (4) Stop bit(s), corresponding to a binary 1.

5.2.2.2 The time between the stop bit of one character and the start bit of the next character may be of any duration. The data interchange circuit is in the marking condition between characters.

5.2.2.3 Even parity corresponds to a parity bit chosen in such a way that there are an even number of ONE bits in the sequence of data bits and parity bit. Odd parity corresponds to an odd number of ONE bits when formed in the same way.

5.2.2.4 All devices must be capable of sending and receiving characters consisting of one start bit, eight data bits, no parity bit, and one stop bit.

5.2.2.5 The default character structure consists of one start bit, eight data bits, no parity bit, and one stop bit. Eight-data-bit character sets are allowed but not specified by this standard. Other character structures can be used for specialized applications, for example, seven data bits, odd, even, mark or space parity, or two stop bits.

5.2.2.6 The character bit sequencing, structure, and parity sense definitions conform to ANSI standards X3.15-1976 and X3.16-1976.<sup>12,13</sup>


### 5.2.3 Speed

5.2.3.1 The data transmission rate for instruments shall be at least one of these baud rates: 1200, 2400, 4800, or 9600 baud. The preferred rate is 9600 baud and should be the default setting of the instrument when more than one baud rate is available. The computer system must have the capability for all four baud rates.


5.2.3.2 Devices may optionally have the capability for other baud rates, such as 300, 19 200, and 38 400 baud for use in specialized applications.

### 5.2.4 Interface Connections

5.2.4.1 The conforming connection specified here defines the point of interconnection between the domain of the instrument and domain of the computer system. (See Figures 1 and 2.) Within the domain of either device, any appropriate connection system may be used, preferably with suitable cable locking hardware.



**Figure 1. Connector Strategy for Instrument Computer Connection—Cable Mounted**



**Figure 2. Connector Strategy for Instrument Computer Connections—Chassis Mounted**

5.2.4.2 The conforming connection utilizes a 25-position connector. The connector contact assignments are listed in Table 1. Connector contacts not listed are unused. The connector contact assignments conform to the EIA-232-D-1986 standard for the circuits that are used.<sup>11</sup>

**Table 1. Connector Contact Assignments for RS 232 D Serial Binary Data Exchange**

Contact No.	EIA Circuit	Description	Direction	
			Instrument	Computer
1	...	Shield	...	No Connection
2	BA	Transmitted Data	Output	Input
3	BB	Received Data	Input	Output
7	AB	Signal Ground	...	...

- 5.2.4.3 Contact 1 is the shield connection; it connects to the instrument's (the DTE) frame. The shield connection is left open at the computer (the DCE) to avoid ground loops. There will be no connections on any other pins. All other pins will be open circuits.

## 5.3 Mechanical Characteristics

### 5.3.1 Connector

- 5.3.1.1 The conforming connector associated with the instrument is a commercial type DB-25P (subminiature D male) style connector. The conforming connector associated with the computer is a commercial type DB-25S (subminiature D female) style connector. The connector dimensions must correspond to those given in the EIA-232-D-1986 standard.<sup>11</sup>
- 5.3.1.2 When the conforming connector of the instrument is cable mounted, it shall be configured with a locking device, such as No. 4-40 or M-3 thread female screw locking hardware. When the conforming connector of the computer is cable mounted, it shall be configured with a locking device, such as No. 4-40 or M-3 thread male screw locking hardware. (See Figure 1.)
- 5.3.1.3 When the conforming connector of either device is chassis mounted, it shall be configured with devices such as No. 4-40 or M-3 thread female screw locking hardware. The mating cable connector shall use devices such as No. 4-40 or M-3 thread male screw locking hardware. (See Figure 2.)
- 5.3.1.4 When the conforming connector of the instrument is cable mounted and the conforming connector of the computer is chassis mounted, then a change in the cable-mounted locking hardware is necessary.

### 5.3.2 Cable

Any extension cables to connect the instrument to the computer require a female connector on one end to mate with the instrument and a male connector on the other end to mate with the computer. Detailed requirements of an interconnecting cable are undefined, but good engineering practice should be followed in selecting the cable and connectors. Shielded cable and connectors may be necessary to suppress electromagnetic interface (EMI). Low-capacitance cable may be necessary for long cable lengths or the higher data rates. Appropriate connector-locking hardware should be used at the conforming connectors.

## 6 Data Link Layer for Serial Binary Data Exchange

### 6.1 Overview

The data link layer has procedures for link connection and release; delimiting and synchronism; sequential control; error detection; and error recovery.

- 6.1.1 Link connection and release establish which system sends and which system receives information. Delimiting and synchronism provide for framing of the data and recognition of frames. Sequence control maintains the sequential order of information across the connection. Error detection senses transmission or format errors. Error recovery attempts to recover from detected errors by retransmitting defective frames or returning the link to a neutral state from otherwise unrecoverable errors.
- 6.1.2 The data link layer uses a character-oriented protocol to send messages between directly connected systems. (See ANSI X3.4-1986. Also, see Appendix B for the coding of the ASCII characters.<sup>14</sup>) Some restrictions are placed on the characters that can appear in the message content.

**6.1.3** The data link mode of operation is one-way transfer of information with alternate supervision. Information flows in one direction at a time. Replies occur after information is sent, never at the same time. It is a simplex stop-and-wait protocol.

**6.1.4** At times, the two systems are actively operating to transfer information. The remainder of the time the data link is in a neutral state. See the state diagram in Appendix A.

**6.1.5** There are three distinct phases in transferring information between instrument and computer system. In each phase, one system directs the operation and is responsible for continuity of the communication. The three phases ensure the actions of the sender and receiver are coordinated. The three phases are establishment, transfer, and termination.

## **6.2 Establishment Phase (Link Connection)**

**6.2.1** The establishment phase determines the direction of information flow and prepares the receiver to accept information.

**6.2.2** The sender notifies the receiver that information is available. The receiver responds that it is prepared to receive before information is transmitted.

**6.2.3** A system that does not have information to send normally monitors the data link to detect the establishment phase. It acts as a receiver, waiting for the other system.

**6.2.4** The system with information available initiates the establishment phase. After the sender determines the data link is in a neutral state, it transmits the <ENQ> transmission control character to the intended receiver. The sender will ignore all responses other than <ACK>, <NAK>, or <ENQ>.

**6.2.5** Upon receiving the <ENQ>, the receiver prepares to receive information. All other characters are ignored. It replies with the <ACK> transmission control character signifying it is ready. With this sequence of events, the establishment phase ends and the transfer phase begins.

**6.2.6** A receiver that cannot immediately receive information replies with the <NAK> transmission control character. Upon receiving <NAK>, the sender must wait at least ten seconds before transmitting another <ENQ>.

**6.2.7** Systems not having the ability to receive information always respond to an <ENQ> by replying with a <NAK>. Systems not having the ability to send information never transmit an <ENQ>.

### **6.2.7.1 Contention**

Should both systems simultaneously transmit an <ENQ>, the data link is in contention. The instrument system has priority to transmit information when contention occurs. Contention is resolved as follows:

(1) Upon receiving a reply of <ENQ> to its transmitted <ENQ>, the computer system must stop trying to transmit; it must prepare to receive. When the next <ENQ> is received, it replies with an <ACK> or <NAK>, depending on its readiness to receive.

(2) Upon receiving a reply of <ENQ> to its transmitted <ENQ>, the instrument must wait at least one second before sending another <ENQ>.



## 6.3 Transfer Phase

During the transfer phase, the sender transmits messages to the receiver. The transfer phase continues until all messages are sent.

### 6.3.1 Frames

Messages are sent in frames; each frame contains a maximum of 64 000 characters (including frame overhead). Messages longer than 64 000 characters are divided between two or more frames.

6.3.1.1 Multiple messages are never combined in a single frame. Every message must begin in a new frame.

6.3.1.2 A frame is one of two types: an intermediate frame or an end frame. Intermediate frames terminate with the characters <ETB>, checksum, <CR>, and <LF>. End frames terminate with the characters <ETX>, checksum, <CR>, and <LF>. A message containing 64 000 characters or less is sent in a single end frame. Longer messages are sent in intermediate frames, with the last part of the message sent in an end frame. The frame structure is illustrated as follows:

<STX> FN text <ETB> C1 C2 <CR> <LF> ← intermediate frame  
 <STX> FN text <ETX> C1 C2 <CR> <LF> ← end frame

where:

<STX>	=	Start of Text transmission control character
FN	=	single digit Frame Number <b>0</b> to <b>7</b>
text	=	Data Content of Message
<ETB>	=	End of Transmission Block transmission control character
<ETX>	=	End of Text transmission control character
C1	=	most significant character of checksum <b>0</b> to <b>9</b> and <b>A</b> to <b>F</b>
C2	=	least significant character of checksum <b>0</b> to <b>9</b> and <b>A</b> to <b>F</b>
<CR>	=	Carriage Return ASCII character
<LF>	=	Line Feed ASCII character

### 6.3.2 Frame Number

The frame number permits the receiver to distinguish between new and retransmitted frames. It is a single digit sent immediately after the <STX> character.

6.3.2.1 The frame number is an ASCII digit ranging from **0** to **7**. The frame number begins at **1** with the first frame of the transfer phase. The frame number is incremented by one for every new frame transmitted. After **7**, the frame number rolls over to **0** and continues in this fashion.

### 6.3.3 Checksum

The checksum permits the receiver to detect a defective frame. The checksum is encoded as two characters, which are sent after the <ETB> or <ETX> character. The checksum is computed by adding the binary values of the characters, keeping the least significant eight bits of the result.

- 6.3.3.1 The checksum is initialized to zero with the <STX> character. The first character used in computing the checksum is the frame number. Each character in the message text is added to the checksum (modulo 256). The computation for the checksum does not include <STX>, the checksum characters, or the trailing <CR> and <LF>.
- 6.3.3.2 The checksum is an integer represented by eight bits; it can be considered as two groups of four bits. The groups of four bits are converted to the ASCII characters of the hexadecimal representation. The two ASCII characters are transmitted as the checksum, with the most significant character first.
- 6.3.3.3 For example, a checksum of 122 can be represented as 01111010 in binary or **7A** in hexadecimal. The checksum is transmitted as the ASCII character **7** followed by the character **A**.

#### 6.3.4 Acknowledgments

After a frame is sent, the sender stops transmitting until a reply is received.

- 6.3.4.1 The receiver replies to each frame. When it is ready to receive the next frame, it transmits one of three replies to acknowledge the last frame. This reply must be transmitted within the timeout period specified in Section 6.5.2.
- 6.3.4.2 A reply of <ACK> signifies the last frame was received successfully, and the receiver is prepared to receive another frame. The sender must increment the frame number and either send a new frame or terminate.
- 6.3.4.3 A reply of <NAK> signifies the last frame was not successfully received, and the receiver is prepared to receive the frame again.
- 6.3.4.4 A reply of <EOT> signifies the last frame was received successfully; the receiver is prepared to receive another frame but requests that the sender stop transmitting. (See the following section on receiver interrupts.)

#### 6.3.5 Receiver Interrupts

The receiver interrupt is a means for the receiver to request the sender to stop transmitting messages as soon as possible.

- 6.3.5.1 During the transfer phase, if the receiver responds to a frame with an <EOT> in place of the usual <ACK>, the sender must interpret this reply as a receiver interrupt request. The <EOT> is a positive acknowledgment of the end frame, signifies the receiver is prepared to receive next frame, and is a request to the sender to stop transmitting.
- 6.3.5.2 The sender does not have to stop transmitting after receiving the receiver interrupt request. If the sender chooses to ignore the <EOT>, the receiver must re-request the interrupt for the request to remain valid.
- 6.3.5.3 If the sender chooses to honor the receiver interrupt request, it must first enter the termination phase to return the data link to the neutral state. This gives the receiver an opportunity to enter the establishment phase and become the sender. The original sender must not enter the establishment phase for at least 15 seconds or until the receiver has sent a message and returned the data link to the neutral state.

## 6.4 Termination Phase (Link Release)

The termination phase returns the data link to the clear or neutral state. The sender notifies the receiver that all messages have been sent.

**6.4.1** The sender transmits the <EOT> transmission control character and then regards the data link to be in a neutral state. Upon receiving <EOT>, the receiver also regards the data link to be in the neutral state.

## 6.5 Error Recovery

Methods are described that enable both sender and receiver to recover, in an orderly way, from errors in data transmission.

### 6.5.1 Defective Frames

A receiver checks every frame to guarantee it is valid. A reply of <NAK> is transmitted for invalid frames. Upon receiving the <NAK>, the sender retransmits the last frame with the same frame number. In this way, transmission errors are detected and automatically corrected.

**6.5.1.1** Any characters occurring before the <STX> or <EOT> or after the end of the block character (the <ETB> or <ETX>) are ignored by the receiver when checking the frame. A frame should be rejected because:

- (1) Any character errors are detected (parity error, framing error, etc.);
- (2) The frame checksum does not match the checksum computed on the received frame; or
- (3) The frame number is not the same as the last accepted frame or one number higher (modulo 8).

**6.5.1.2** Upon receiving a <NAK> or any character except an <ACK> or <EOT> (a <NAK> condition), the sender increments a retransmit counter and retransmits the frame. If this counter shows a single frame was sent and not accepted six times, the sender must abort this message by proceeding to the termination phase. An abort should be extremely rare, but it provides a mechanism to escape from a condition where the transfer phase cannot continue.

### 6.5.2 Timeouts

The sender and receiver both use timers to detect loss of coordination between them. The timers provide a method for recovery if the communication line or the other device fails to respond.

**6.5.2.1** During the establishment phase, the sender sets a timer when transmitting the <ENQ>. If a reply of an <ACK>, <NAK>, or <ENQ> is not received within 15 seconds, a timeout occurs. After a timeout, the sender enters the termination phase.

**6.5.2.2** During the establishment phase, if the computer (as receiver) detects contention, it sets a timer. If an <ENQ> is not received within 20 seconds, a timeout occurs. After a timeout, the receiver regards the line to be in the neutral state.

**6.5.2.3** During the transfer phase, the sender sets a timer when transmitting the last character of a frame. If a reply is not received within 15 seconds, a timeout occurs. After a timeout, the sender aborts the message transfer by proceeding to the termination phase. As with excessive retransmissions of defective frames, the message must be remembered so it can be completely repeated.

- 6.5.2.4 During the transfer phase, the receiver sets a timer when first entering the transfer phase or when replying to a frame. If a frame or <EOT> is not received within 30 seconds, a timeout occurs. After a timeout, the receiver discards the last incomplete message and regards the line to be in the neutral state.
- 6.5.2.5 A receiver must reply to a frame within 15 seconds or the sender will timeout. A receiver can delay its reply for up to 15 seconds to process the frame or to otherwise go into “busy” mode. Longer delays cause the sender to abort the message.
- 6.5.2.6 Receivers that cannot process messages fast enough to keep up with a sender may cause message buffer overflows in the sender. A sender can normally store at least one complete message. Storage space for more than one outgoing message is desirable but optional.

## 6.6 Restricted Message Characters

The data link protocol is designed for sending character-based message text. Restrictions are placed on which characters may appear in the message text. The restrictions make it simpler for senders and receivers to recognize replies and frame delimiters. Additional characters are restricted to avoid interfering with software controls for devices such as multiplexers.

- 6.6.1 A <LF> character is not permitted to appear in the message text; it can appear only as the last character of a frame.
- 6.6.2 None of the ten transmission control characters, the <LF> format effector control character, or four device control characters may appear in message text. The restricted characters are: <SOH>, <STX>, <ETX>, <EOT>, <ENQ>, <ACK>, <DLE>, <NAK>, <SYN>, <ETB>, <LF>, <DC1>, <DC2>, <DC3>, and <DC4>.

## 7 Physical Layer for TCP/IP Data Exchange

### 7.1 Overview

The mechanical and electrical connection for TCP/IP data transmission between instrument and computer system is described in this section. The topology is not a point-to-point, direct connection between two devices, but utilizes a central hub for device connection.

### 7.2 Electrical Characteristics

The voltage and impedance levels for the generator and receiver circuits are as specified in the IEEE 802.3 standard.<sup>15</sup>

#### 7.2.1 Signal Levels

- 7.2.1.1 The signal levels conform to the IEEE 802.3 standard.<sup>15</sup>

#### 7.2.2 Character Structure

- 7.2.2.1 The method of data transmission is TCP/IP. The order of the bits in a character is:

- (1) One start bit, corresponding to a binary 0;
- (2) The data bits of the character, least significant bit transmitted first;
- (3) Parity bit; and
- (4) Stop bit(s), corresponding to a binary 1.


- 7.2.2.2 The time between the stop bit of one character and the start bit of the next character may be of any duration. The data interchange circuit is in the marking condition between characters.
- 7.2.2.3 Even parity corresponds to a parity bit chosen in such a way that there are an even number of ONE bits in the sequence of data bits and parity bit. Odd parity corresponds to an odd number of ONE bits when formed in the same way.
- 7.2.2.4 All devices must be capable of sending and receiving characters consisting of one start bit, eight data bits, no parity bit, and one stop bit.
- 7.2.2.5 The default character structure consists of one start bit, eight data bits, no parity bit, and one stop bit. Eight-data-bit character sets are allowed but not specified by this standard. Other character structures can be used for specialized applications, for example, seven data bits, odd, even, mark or space parity, or two stop bits.
- 7.2.2.6 The character bit sequencing, structure, and parity sense definitions conform to ANSI standards X3.15-1976 and X3.16-1976.<sup>12,13</sup>

### 7.2.3 Speed

- 7.2.3.1 The data transmission rate for instruments shall conform to IEEE 802.3<sup>15</sup> and operate at least 10 MB/second. A computer system using TCP/IP must have the capability to conform to a minimum speed of 10 MB/second.


### 7.2.4 Interface Connections

- 7.2.4.1 The conforming connection specified here defines the point of interconnection between the domain of the instrument and domain of the computer system. (See Figures 3 and 4.) Within the domain of the devices, an RJ45 connection system is used.



**NOTE:** Looking at connector end with cable running away from you.

**Figure 3. Plug**



**NOTE:** Looking at cavity in the wall.

**Figure 4. Jack**

- 7.2.4.2 The conforming connection utilizes an eight-line connector. The connector contact assignments are listed in Table 2. The connector contact assignments conform to the ANSI EIA/TIA 568B standard<sup>16</sup> (also called the AT&T specification).

**Table 2. Connector Contact Assignments**

Pin	Color
1	white/orange
2	orange
3	white/green
4	blue
5	white/blue
6	green
7	white/brown
8	brown

### 7.3 Mechanical Characteristics

#### 7.3.1 Connector

7.3.1.1 The conforming connector associated with the instrument is a commercial type RJ-45F style connector. The conforming connector associated with the computer is a commercial type RJ-45F style connector. The connector dimensions must correspond to those given in the ANSI EIA/TIA 568B standard.<sup>16</sup>

#### 7.3.2 Cable

Category 5 cable as defined by ANSI EIA/TIA 568B<sup>16</sup> is preferred for all connections. In general, no extension cables should be required to connect the instrument to the computer network. Detailed requirements of an interconnecting cable are undefined, but good engineering practice should be followed in selecting the cable and connectors. Low-capacitance cable and shielded connectors may be necessary to suppress electromagnetic interface (EMI). Appropriate connector-locking hardware should be used at the conforming connectors.

## 8 Data Link Layer for TCP/IP Data Exchange

### 8.1 Overview

The data link layer has procedures for link connection and release; delimiting and synchronism; sequential control; error detection; and error recovery.

**8.1.1** Link connection and release establish which system sends and which system receives information. Delimiting and synchronism provide for framing of the data and recognition of frames. Sequence control maintains the sequential order of information across the connection. Error detection senses transmission or format errors. Error recovery attempts to recover from detected errors by retransmitting defective messages or returning the link to a neutral state from otherwise unrecoverable errors.

**8.1.2** The data link layer uses a character-oriented protocol to send messages between directly connected systems. (See ANSI X3.4-1986. Also, see Appendix B for the coding of the ASCII characters.<sup>14</sup>) Some restrictions are placed on the characters that can appear in the message content.

**8.1.3** The data link mode of operation is one-way transfer of information with alternate supervision. Information flows in one direction at a time. Replies occur after information is sent, never at the same time. It is a simplex stop-and-wait protocol.

**8.1.4** At times, the two systems are actively operating to transfer information. The remainder of the time the data link is in a neutral state. See the state diagram in Appendix A.

**8.1.5** There are three distinct phases in transferring information between instrument and computer systems. In each phase, one system directs the operation and is responsible for continuity of the communication. The three phases ensure the actions of the sender and receiver are coordinated. The three phases are establishment, transfer, and termination.

## **8.2 Establishment Phase (Link Connection)**

**8.2.1** The establishment phase determines the direction of information flow and prepares the receiver to accept information.

**8.2.1.1** The standard framework for TCP/IP network implementation is generally referred to as a “socket” implementation. The normal implementation of a socket-based communication protocol is asymmetric. One end is considered a “server” providing service to the other end, which is considered a “client.” Under this client/server paradigm, the “computer system” is defined as the “server,” and the “instrument” is defined as the “client.” Both ends must initially create sockets and bind their local Internet address to their socket. The computer system, as the server, offers a socket for connection. The instrument (client) must be independently configurable for both socket and IP address (to facilitate use of server-defined socket number and IP address), so it can request a connection to the offered socket.

**8.2.2** The sender notifies the receiver that information is available. The receiver responds that it is prepared to receive before information is transmitted.

**8.2.3** A system that does not have information to send normally monitors the data link to detect the establishment phase. It acts as a receiver, waiting for the other system.

**8.2.4** The system with information available initiates the establishment phase. After the sender determines the data link is in a neutral state, it transmits the <ENQ> transmission control character to the intended receiver. The sender will ignore all responses other than <ACK>, <NAK>, or <ENQ>.

**8.2.5** Upon receiving the <ENQ>, the receiver prepares to receive information. All other characters are ignored. It replies with the <ACK> transmission control character signifying it is ready. With this sequence of events, the establishment phase ends and the transfer phase begins.

**8.2.6** A receiver that cannot immediately receive information replies with the <NAK> transmission control character. Upon receiving <NAK>, the sender must wait at least ten seconds before transmitting another <ENQ>.

**8.2.7** Systems not having the ability to receive information always respond to an <ENQ> by replying with a <NAK>. Systems not having the ability to send information never transmit an <ENQ>.

### **8.2.7.1 Contention**

Should both systems simultaneously transmit an <ENQ>, the data link is in contention. The instrument system has priority to transmit information when contention occurs. Contention is resolved as follows:

- (1) Upon receiving a reply of <ENQ> to its transmitted <ENQ>, the computer system must stop trying to transmit; it must prepare to receive. When the next <ENQ> is received, it replies with an <ACK> or <NAK>, depending on its readiness to receive.



- (2) Upon receiving a reply of <ENQ> to its transmitted <ENQ>, the instrument must wait at least one second before sending another <ENQ>.

### 8.3 Transfer Phase

During the transfer phase, the sender transmits messages to the receiver. The transfer phase continues until all messages are sent.

#### 8.3.1 Frames

Messages are sent in frames; each frame contains a maximum of 64 000 characters (including frame overhead). Messages longer than 63 993 characters are divided between two or more frames.

8.3.1.1 Multiple messages are never combined in a single frame. Every message must begin in a new frame.

8.3.1.2 A frame is one of two types: an intermediate frame or an end frame. Intermediate frames terminate with the characters <ETB>, checksum, <CR>, and <LF>. End frames terminate with the characters <ETX>, checksum, <CR>, and <LF>. A message containing 63 993 characters or less is sent in a single end frame. Longer messages are sent in intermediate frames, with the last part of the message sent in an end frame. The frame structure is illustrated as follows:

<STX> FN text <ETB> C1 C2 <CR> <LF> ← intermediate frame

<STX> FN text <ETX> C1 C2 <CR> <LF> ← end frame

where:

<STX>	=	Start of Text transmission control character
FN	=	single digit Frame Number <b>0</b> to <b>7</b>
text	=	Data Content of Message
<ETB>	=	End of Transmission Block transmission control character
<ETX>	=	End of Text transmission control character
C1	=	most significant character of checksum <b>0</b> to <b>9</b> and <b>A</b> to <b>F</b>
C2	=	least significant character of checksum <b>0</b> to <b>9</b> and <b>A</b> to <b>F</b>
<CR>	=	Carriage Return ASCII character
<LF>	=	Line Feed ASCII character

#### 8.3.2 Frame Number

The frame number permits the receiver to distinguish between new and retransmitted frames. It is a single digit sent immediately after the <STX> character.

8.3.2.1 The frame number is an ASCII digit ranging from **0** to **7**. The frame number begins at **1** with the first frame of the transfer phase. The frame number is incremented by one for every new frame transmitted. After **7**, the frame number rolls over to **0** and continues in this fashion.



### 8.3.3 Checksum

The checksum permits the receiver to detect a defective frame. The checksum is encoded as two characters, which are sent after the <ETB> or <ETX> character. The checksum is computed by adding the binary values of the characters, keeping the least significant eight bits of the result.

8.3.3.1 The checksum is initialized to zero with the <STX> character. The first character used in computing the checksum is the frame number. Each character in the message text is added to the checksum (modulo 256). The computation for the checksum does not include <STX>, the checksum characters, or the trailing <CR> and <LF>.

8.3.3.2 The checksum is an integer represented by eight bits; it can be considered as two groups of four bits. The groups of four bits are converted to the ASCII characters of the hexadecimal representation. The two ASCII characters are transmitted as the checksum, with the most significant character first.

8.3.3.3 For example, a checksum of 122 can be represented as 01111010 in binary or **7A** in hexadecimal. The checksum is transmitted as the ASCII character **7** followed by the character **A**.

### 8.3.4 Acknowledgments

After a frame is sent, the sender stops transmitting until a reply is received.

8.3.4.1 The receiver replies to each frame. When it is ready to receive the next frame, it transmits one of three replies to acknowledge the last frame. This reply must be transmitted within the timeout period specified in Section 8.5.2.

8.3.4.2 A reply of <ACK> signifies the last frame was received successfully and the receiver is prepared to receive another frame. The sender must increment the frame number and either send a new frame or terminate.

8.3.4.3 A reply of <NAK> signifies the last frame was not successfully received and the receiver is prepared to receive the frame again.

8.3.4.4 A reply of <EOT> signifies the last frame was received successfully, the receiver is prepared to receive another frame, but requests that the sender stop transmitting. (See the following section on receiver interrupts.)

### 8.3.5 Receiver Interrupts

The receiver interrupt is a means for the receiver to request the sender to stop transmitting messages as soon as possible.

8.3.5.1 During the transfer phase, if the receiver responds to a frame with an <EOT> in place of the usual <ACK>, the sender must interpret this reply as a receiver-interrupt request. The <EOT> is a positive acknowledgment of the end frame, signifies the receiver is prepared to receive next frame, and is a request to the sender to stop transmitting.

8.3.5.2 The sender does not have to stop transmitting after receiving the receiver-interrupt request. If the sender chooses to ignore the <EOT>, the receiver must re-request the interrupt for the request to remain valid.

8.3.5.3 If the sender chooses to honor the receiver-interrupt request, it must first enter the termination phase to return the data link to the neutral state. This gives the receiver an opportunity to enter the establishment phase and become the sender. The original sender must not enter the establishment phase for at least 15 seconds or until the receiver has sent a message and returned the data link to the neutral state.

## 8.4 Termination Phase (Link Release)

The termination phase returns the data link to the clear or neutral state. The sender notifies the receiver that all messages have been sent.

8.4.1 The sender transmits the <EOT> transmission control character and then regards the data link to be in a neutral state. Upon receiving <EOT>, the receiver also regards the data link to be in the neutral state.

## 8.5 Error Recovery

Methods are described that enable both sender and receiver to recover, in an orderly way, from errors in data transmission.

### 8.5.1 Defective Frames

A receiver checks every frame to guarantee it is valid. A reply of <NAK> is transmitted for invalid frames. Upon receiving the <NAK>, the sender retransmits the last frame with the same frame number. In this way, transmission errors are detected and automatically corrected.

8.5.1.1 Any characters occurring before the <STX> or <EOT> or after the end of the block character (the <ETB> or <ETX>) are ignored by the receiver when checking the frame. A frame should be rejected because:

- (1) Any character errors are detected (parity error, framing error, etc.);
- (2) The frame checksum does not match the checksum computed on the received frame; or
- (3) The frame number is not the same as the last accepted frame or one number higher (modulo 8).

8.5.1.2 Upon receiving a <NAK> or any character except an <ACK> or <EOT> (a <NAK> condition), the sender increments a retransmit counter and retransmits the frame. If this counter shows a single frame was sent and not accepted six times, the sender must abort this message by proceeding to the termination phase. An abort should be extremely rare, but it provides a mechanism to escape from a condition where the transfer phase cannot continue.

### 8.5.2 Timeout

The sender and receiver both use timers to detect loss of coordination between them. The timers provide a method for recovery if the communication line or the other device fails to respond.

8.5.2.1 During the establishment phase, the sender sets a timer when transmitting the <ENQ>. If a reply of an <ACK>, <NAK>, or <ENQ> is not received within 15 seconds, a timeout occurs. After a timeout, the sender enters the termination phase.

8.5.2.2 During the establishment phase, if the computer (as receiver) detects contention, it sets a timer. If an <ENQ> is not received within 20 seconds, a timeout occurs. After a timeout, the receiver regards the line to be in the neutral state.

- 8.5.2.3 During the transfer phase, the sender sets a timer when transmitting the last character of a frame. If a reply is not received within 15 seconds, a timeout occurs. After a timeout, the sender aborts the message transfer by proceeding to the termination phase. As with excessive retransmissions of defective frames, the message must be remembered so it can be completely repeated.
- 8.5.2.4 During the transfer phase, the receiver sets a timer when first entering the transfer phase or when replying to a frame. If a frame or <EOT> is not received within 30 seconds, a timeout occurs. After a timeout, the receiver discards the last incomplete message and regards the line to be in the neutral state.
- 8.5.2.5 A receiver must reply to a frame within 15 seconds or the sender will timeout. A receiver can delay its reply for up to 15 seconds to process the frame or to otherwise go into “busy” mode. Longer delays cause the sender to abort the message.
- 8.5.2.6 Receivers that cannot process messages fast enough to keep up with a sender may cause message buffer overflows in the sender. A sender can normally store at least one complete message. Storage space for more than one outgoing message is desirable but optional.

## 8.6 Restricted Message Characters

The data link protocol is designed for sending character-based message text. Restrictions are placed on which characters may appear in the message text. The restrictions make it simpler for senders and receivers to recognize replies and frame delimiters. Additional characters are restricted to avoid interfering with software controls for devices such as multiplexers.

- 8.6.1 A <LF> character is not permitted to appear in the message text; it can appear only as the last character of a frame.
- 8.6.2 None of the ten transmission control characters, the <LF> format effector control character, or four device control characters may appear in message text. The restricted characters are: <SOH>, <STX>, <ETX>, <EOT>, <ENQ>, <ACK>, <DLE>, <NAK>, <SYN>, <ETB>, <LF>, <DC1>, <DC2>, <DC3>, and <DC4>.

## References

- <sup>1</sup> CLSI/NCCLS. *Specification for Transferring Information Between Clinical Laboratory Instruments and Information Systems; Approved Standard—Second Edition*. CLSI/NCCLS document LIS2-A2. Wayne, PA: NCCLS; 2004.
- <sup>2</sup> CLSI/NCCLS. *Standard Guide for Selection of a Clinical Laboratory Information Management System*. CLSI/NCCLS document LIS3-A. Wayne, PA: NCCLS; 2003.
- <sup>3</sup> CLSI/NCCLS. *Standard Guide for Documentation of Clinical Laboratory Computer Systems*. CLSI/NCCLS document LIS4-A. Wayne, PA: NCCLS; 2003.
- <sup>4</sup> CLSI/NCCLS. *Standard Specification for Transferring Clinical Observations Between Independent Computer Systems*. CLSI/NCCLS document LIS5-A. Wayne, PA: NCCLS; 2003.
- <sup>5</sup> CLSI/NCCLS. *Standard Practice for Reporting Reliability of Clinical Laboratory Information Systems*. CLSI/NCCLS document LIS6-A. Wayne, PA: NCCLS; 2003.
- <sup>6</sup> CLSI/NCCLS. *Standard Specification for Use of Bar Codes on Specimen Tubes in the Clinical Laboratory*. CLSI/NCCLS document LIS7-A. Wayne, PA: NCCLS; 2003.
- <sup>7</sup> CLSI/NCCLS. *Standard Guide for Functional Requirements of Clinical Laboratory Information Management Systems*. CLSI/NCCLS document LIS8-A. Wayne, PA: NCCLS; 2003.
- <sup>8</sup> CLSI/NCCLS. *Standard Guide for Coordination of Clinical Laboratory Services within the Electronic Health Record Environment and Networked Architectures*. CLSI/NCCLS document LIS9-A. Wayne, PA: NCCLS; 2003.
- <sup>9</sup> ISO. *Information processing systems—Open systems interconnection—Basic reference model*. ISO 7498-1984(E). Geneva: International Organization for Standardization
- <sup>10</sup> CLSI/NCCLS. *Laboratory Automation: Communications with Automated Clinical Laboratory Systems, Instruments, Devices, and Information Systems; Approved Standard*. CLSI/NCCLS document AUTO3-A. Wayne, PA: NCCLS; 2000.
- <sup>11</sup> EIA. *Interface Between Data Terminal Equipment and Data Circuit-Terminating Equipment Employing Serial Binary Data Interchange*. EIA-232-D-1986. Washington, DC: Electronics Industry Association; 1986.
- <sup>12</sup> ANSI. *American National Standard for Bit Sequencing of the American National Standard Code for Information Interchange in Serial-by-Bit Data Transmission*. X3.15-1976. New York, NY: American National Standards Institute; 1976.
- <sup>13</sup> ANSI. *American National Standard Character Structure and Character Parity Sense for Serial-by-Bit Data Communication in the American National Standard Code for Information Interchange*. X3.16-1976. New York, NY: American National Standards Institute; 1976.
- <sup>14</sup> ANSI. *American National Standard Code for Information Systems—Coded Character Sets—7-Bit American National Standard Code for Information Interchange (7-Bit ASCII)*. X3.4-1986. New York, NY: American National Standards Institute; 1986.
- <sup>15</sup> IEEE. *Information Technology—Telecommunications and Information Exchange Between Systems—Local and Metropolitan Area Networks—Specific Requirements—Part 3: Carrier Sense Multiple Access With Collision Detection (CSMA/CD) Access Method and Physical Layer Specifications*. IEEE-802.3. Piscataway, NJ: Institute of Electrical and Electronics Engineers; 2000.
- <sup>16</sup> ANSI EIA/TIA 568B. *Commercial Building Telecommunications Cabling Standard*. Arlington, VA: Electronic Industries Alliance; 2001.

## Additional Reference

INCITS. *Information Technology – ISO 7-bit coded Character Set for Information Interchange [L2]*. INCITS/ISO/IEC 646:1991[R2004]. Washington, DC: InterNational Committee for Information Technology Standards; 2004.

## Appendix A. Mandatory Information




Figure A1. State Diagram

NOTE 1: “%8” represents modulo 8.

NOTE 2: “=” represents assignment of a value. “Timer: = 15” resets the timer to 15 s as used here.

## B1. Seven-Bit ASCII Code Charts

[illegible]

**Figure B1. Decimal Character Code**

**| hex CHR |**

↑      ↑  
**\_ASCII Character\_**  
**-Hexadecimal Character Code-**

00 NUL	10 DLE	20 SP	30 0	40 @	50 P	60 `	70 p
01 SOH	11 DC1	21 !	31 1	41 A	51 Q	61 a	71 q
02 STX	12 DC2	22 "	32 2	42 B	52 R	62 b	72 r
03 ETX	13 DC3	23 #	33 3	43 C	53 S	63 c	73 s
04 EOT	14 DC4	24 \$	34 4	44 D	54 T	64 d	74 t
05 ENQ	15 NAK	25 %	35 5	45 E	55 U	65 e	75 u
06 ACK	16 SYN	26 &	36 6	46 F	56 V	66 f	76 v
07 BEL	17 ETB	27 '	37 7	47 G	57 W	67 g	77 w
08 BS	18 CAN	28 (	38 8	48 H	58 X	68 h	78 x
09 HT	19 EM	29 )	39 9	49 I	59 Y	69 i	79 y
0A LF	1A SUB	2A *	3A :	4A J	5A Z	6A j	7A z
0B VT	1B ESC	2B +	3B ;	4B K	5B [	6B k	7B {
0C FF	1C FS	2C ,	3C <	4C L	5C \	6C l	7C
0D CR	1D GS	2D -	3D =	4D M	5D ]	6D m	7D }
0E SO	1E RS	2E .	3E >	4E N	5E ^	6E n	7E ~
0F SI	1F US	2F /	3F ?	4F O	5F _	6F o	7F DEL

### Figure B2. Hexadecimal Character Code

**Clinical and Laboratory Standards Institute consensus procedures include an appeals process that is described in detail in Section 8 of the Administrative Procedures. For further information, contact CLSI or visit our website at [www.clsi.org](http://www.clsi.org).**

## Summary of Delegate Comments and Area Committee Responses

LIS01-A2: *Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems; Approved Standard—Second Edition*

### General

1. This is a review document of the RS-232 interface, which is rarely used now. Most laboratory equipment comes as an integrated system with simple interfacing software.
  - **There are still IVD, middleware, and LIS vendors that make use of this specification for interfaces to legacy and new IVD instruments. There are current IVD vendors that still support RS-232 connectivity to their instruments. Often, this native RS-232 interface may be connected by a very short cable to a device server, either single or multiple ports, but the native instrument datastream is RS-232. The specification still supports these devices.**
2. The College of American Pathologists recommends that the document drafters might want to consider adding commentary about wireless communications possibilities, as the document assumes wired connections.
  - **See the response to comment 7 below.**
3. The document describes Ethernet with TCP/IP as an alternative choice for serial communication. The CLIMS shall be TCP/IP-server, the analyzer TCP/IP-client. There are analyzers out in the field that offer ASTM-communication by TCP/IP, but can only play the TCP/IP-server. There is a need for the CLSI standard to permit analyzers as TCP/IP-servers. This behavior is not wrong because analyzers are servers for diagnostic measurements, too. It is also reasonable for the CLIMS to reorganize the analyzer connection after reboot or network shutdown.
  - **There are advantages to having the CLIMS assigned the server role. Logically a client-server approach implies one server and many clients, which is true for the case of one CLIMS and multiple instruments. Assigning the instrument the role of client puts control of the initiation and maintenance of the TCP/IP session in the hands of what is presumably the less powerful system. Also, the configuration of the network environment can be simpler if only one server must be identified on the network (the CLIMS), rather than multiple servers if each analyzer were a server. In addition, having the analyzer assigned as the client provides a level of protection for the analyzer, as it must initiate communication rather than receive an unsolicited connection. Finally, if the analyzer is protected by a component such as a firewall, the component would need to be configured to allow the inbound traffic if the analyzer is the server. Configuration is not necessary if the analyzer initiates the communication as a client. This is addressed in Section 8.2.1.1.**
4. Most of the document describes the connectivity hardware. In fact, this hardware is described in other standards very well. If LIS01-A2 would refer to these hardware standards without reproducing, then it is updated automatically after hardware improvements. For example, there is described a 25-pin DB-25 S/P connector, but all modern analyzers are fitted with the newer 9-pin connectors. There is a need for 9-pin connectors to be permitted by the standard.
  - **The reviewer makes a good point, but one of the goals for this standard was to maintain compatibility with the previous version. See the response to comment 7 below.**
5. The new standard increases the possible frame length from 255 up to 64 000 characters. This is advantageous for the developer because intermediate frames become no longer necessary, but there is a danger. A 2-digit



checksum cannot hedge 64 000 bytes in case of serial communication. In case of TCP/IP communication, the frame length can be unlimited without danger and even the checksum is obsolete. I think 64 000 is a bad compromise between 256 and infinity, and there is a need for a frame length reduction in case of serial communication.

- **See Sections 4.4 and 4.4.1, which address the reviewer's comments concerning large frame sizes and restricting their use to TCP/IP communication.**
6. In case of TCP/IP communication, frame length limitations and checksums are obsolete; but compatibility with serial communication requires ability of processing checksums and intermediate frames, so the checksum shall be an optional parameter.
  - **As noted by the reviewer, the checksum is no longer required for TCP/IP communication. However, the structure of the frame was not modified for TCP/IP communication in order to maintain complete compatibility with legacy RS-232 interfaces. Compatibility and simplification of frame structure argues against allowing for multiple possible frame structures, and support for an optional checksum parameter introduces opportunities for error based on multiple options. Because the calculation of the checksum is required for serial modes, it introduces no new effort and minimal communication overhead to the standard.**
  7. Since the document has been updated to include TCP/IP communication, would it also be possible to include a small section on other communication technologies, such as virtual COM ports over USB, high speed USB, and firewire (IEEE 1394)?
  - **Because of the nature of the typical laboratory environment, where a central LIS is connected to multiple instruments located at various locations, it is expected the majority of connections will be RS-232 or TCP/IP. The use of virtual COM ports, through the use of devices such as terminal servers, can be implemented in a manner that is transparent to this specification. Likewise, a wireless implementation could be transparent depending on the implementation. This revision was intended to clarify the existing standards, and the implementation of other technologies is beyond the scope of this effort.**
  - **For clarification, Section 4.5 has been added: "Although this specification addresses serial communication using a 25-pin connector and TCP/IP communication, it can be adapted to other standard communication approaches such as a 9-pin serial connector, wireless, high speed USB, and firewire. The physical layer specifications for other approaches must be agreed upon outside the context of this standard. Adaptation to 9-pin, USB, etc. is standard, and can be accomplished in a manner transparent to the specification."**

#### Table of Contents

8. The proposed Table of Contents is less detailed than it was in the first edition. An original construction of the table was more convenient for users. Preserve the original format of the table (from the first edition).
- **Although less detailed, the Table of Contents is consistent with other CLSI documents (such as LIS2-A2). The structure of the Table of Contents should be consistent across all CLSI documents, and questions about the acceptability/usability of that structure should take place at a higher level than this document.**

#### Section 3, Terminology

9. It would be helpful for users to see additional terms that were frequently used in a text. Add descriptions for the following terms: stop bit, parity bit, neutral stage, checksum, instrument domain, and computer domain.
- **There are other terms that are not specifically defined in the terminology section, such as TCP/IP. The specification does reference other standards that provide definitions for most terms not listed in the terminology section. The terms used are common to the communications industry and could reasonably be expected to be understood by a reader/user of this standard, or can be understood by reviewing the referenced standards.**



Section 5.2.3.1

10. The data transmission rate for instruments shall be at least one of these baud rates: 1200, 2400, 4800, or 9600 baud. The preferred rate is 9600 baud, and should be the default setting of the instrument when more than one baud rate is available. The computer system must have the capability for all four baud rates.

Is it still necessary that when faster baud rates are available (ie, 57 600) that the default setting of the instrument be 9600?

Is it still true that instruments must support 1200, 2400, 4800, or 9600? And if so, is this meant to ensure compatibility with older computer systems?

Note: Since the document has been updated to include TCP/IP data exchange, the TCP/IP characteristics specify speed in terms of a “minimum” such as “at least 10MB/sec” vs stating that systems “must support 10MB/sec.”

- **The reviewer makes a good point, but legacy compatibility is an important component of the standard. In order to maintain compatibility with existing implementations, it is best to leave the default value at 9600. Note that the wording does not preclude use of faster transmission rates when possible, but rather requires compatibility to slower computer systems to prevent incompatibilities.**
- **“Minimum” was used so as not to exclude other speeds such as 100MB/sec.**

Section 5.2.4.2

11. The conforming connection utilizes a 25-position connector. The connector contact assignments are listed in Table 1. Connector contacts not listed are unused. The connector contact assignments conform to the EIA-232-D-1986 standard for the circuits that are used.

Is it still a requirement that the conforming interface connection on laboratory instruments utilize a 25-position connector? Considering only four pins are used, and the industry focus on making instruments more compact, a 9-pin connection or possibly a USB type connection would be more appropriate depending on the device. In many cases where serial communication is required, more and more users are running serial-to-Ethernet convertors, or serial-to-USB convertors to attach to systems/networks.

- **See the response to comment 7 above.**

Sections 6.3.1.2 and 8.3.1.2

12. These sections include the same table of frame codes. Include ALL frame code in one comprehensive table that should be added to one of appendixes; exclude tables of frame codes from the body of the standard.
- **The frame layout was included in the original ASTM standard body, as well as the previous CLSI version. The frame layout is included where it might be most immediately referenced, making it easier for a reader to reference the frame definition for a specific implementation (RS-232, TCP/IP). Having to go to an appendix and then back to the body of the document detracts from usability.**

Section 8.2.1.1

13. It might be preferred that the instrument is designated as the server and the computer system is identified as the client. The rationale is that if one has to reboot the computer system for some reason, and it is designated as the “server,” all instruments that are connected to it will need to reconnect to that computer system. If the instruments are “servers,” they are not affected when the computer system is rebooted.

- **See the discussion and response to comment 3 above.**

## The Quality Management System Approach

Clinical and Laboratory Standards Institute (CLSI) subscribes to a quality management system approach in the development of standards and guidelines, which facilitates project management; defines a document structure via a template; and provides a process to identify needed documents. The approach is based on the model presented in the most current edition of CLSI/NCCLS document HS1—*A Quality Management System Model for Health Care*. The quality management system approach applies a core set of “quality system essentials” (QSEs), basic to any organization, to all operations in any health care service’s path of workflow (ie, operational aspects that define how a particular product or service is provided). The QSEs provide the framework for delivery of any type of product or service, serving as a manager’s guide. The QSEs are:

Documents & Records	Equipment	Information Management	Process Improvement
Organization	Purchasing & Inventory	Occurrence Management	Customer Service
Personnel	Process Control	Assessments—External & Internal	Facilities & Safety

LIS01-A2 addresses the QSEs indicated by an “X.” For a description of the other documents listed in the grid, please refer to the Related CLSI Reference Materials section on the following page.

Documents & Records	Organization	Personnel	Equipment	Purchasing & Inventory	Process Control	Information Management	Occurrence Management	Assessments—External & Internal	Process Improvement	Customer Service	Facilities & Safety
					X AUTO3 LIS2	X AUTO3 LIS2					

Adapted from CLSI/NCCLS document HS1—*A Quality Management System Model for Health Care*.

### Path of Workflow

A path of workflow is the description of the necessary steps to deliver the particular product or service that the organization or entity provides. For example, CLSI/NCCLS document GP26—*Application of a Quality Management System Model for Laboratory Services* defines a clinical laboratory path of workflow, which consists of three sequential processes: preexamination, examination, and postexamination. All clinical laboratories follow these processes to deliver the laboratory’s services, namely quality laboratory information.

LIS01-A2 addresses the clinical laboratory path of workflow steps indicated by an “X.” For a description of the other documents listed in the grid, please refer to the Related CLSI Reference Materials section on the following page.

Preexamination				Examination			Postexamination	
Examination ordering	Sample collection	Sample transport	Sample receipt/processing	Examination	Results review and follow-up	Interpretation	Results reporting and archiving	Sample management
							X AUTO3 LIS2	

Adapted from CLSI/NCCLS document HS1—*A Quality Management System Model for Health Care*.

## Related CLSI Reference Materials\*

- AUTO3-A**      **Laboratory Automation: Communications with Automated Clinical Laboratory Systems, Instruments, Devices, and Information Systems; Approved Standard (2000).** This document provides standards to facilitate accurate and timely electronic exchange of data and information between the automated laboratory elements.
- LIS2-A2**      **Specification for Transferring Information Between Clinical Laboratory Instruments and Information Systems; Approved Standard—Second Edition (2004).** This document covers the two-way digital transmission of remote requests and results between clinical laboratory instruments and information systems.

---

\* Proposed-level documents are being advanced through the Clinical and Laboratory Standards Institute consensus process; therefore, readers should refer to the most current editions.

# Active Membership (as of 1 April 2008)

## Sustaining Members

Abbott  
American Association for Clinical Chemistry  
AstraZeneca Pharmaceuticals  
Bayer Corporation  
BD  
Beckman Coulter, Inc.  
bioMérieux, Inc.  
CLMA  
College of American Pathologists  
GlaxoSmithKline  
Ortho-Clinical Diagnostics, Inc.  
Pfizer Inc  
Roche Diagnostics, Inc.

## Professional Members

American Academy of Family Physicians  
American Association for Clinical Chemistry  
American Association for Laboratory Accreditation  
American Association for Respiratory Care  
American Chemical Society  
American College of Medical Genetics  
American Medical Technologists  
American Society for Clinical Laboratory Science  
American Society for Microbiology  
American Type Culture Collection  
ASCP  
Associazione Microbiologi Clinici Italiani (AMCLI)  
Canadian Society for Medical Laboratory Science  
COLA  
College of American Pathologists  
College of Medical Laboratory Technologists of Ontario  
College of Physicians and Surgeons of Saskatchewan  
Elkin Simson Consulting Services  
ESCMID  
Family Health International  
Hong Kong Accreditation Service  
Innovation and Technology Commission  
International Federation of Biomedical Laboratory Science  
International Federation of Clinical Chemistry  
Italian Society of Clinical Biochemistry and Clinical Molecular Biology  
JCCLS  
The Joint Commission  
National Society for Histotechnology, Inc.  
Ontario Medical Association Quality Management Program-Laboratory Service  
RCPA Quality Assurance Programs  
PTY Limited  
Serbian Society of Microbiology  
SIMEI  
Sociedad Espanola de Bioquímica Clínica y Patología Molecular  
Sociedade Brasileira de Analises Clinicas  
Sociedade Brasileira de Patologia Clínica  
Turkish Society of Microbiology  
Washington G2 Reports  
World Health Organization

## Government Members

Armed Forces Institute of Pathology  
Association of Public Health Laboratories  
BC Centre for Disease Control  
Centers for Disease Control and Prevention  
Centers for Disease Control and Prevention – Namibia  
Centers for Disease Control and Prevention – Nigeria  
Centers for Disease Control and Prevention – Tanzania  
Centers for Medicare & Medicaid Services  
Centers for Medicare & Medicaid Services/CLIA Program  
Chinese Committee for Clinical Laboratory Standards

Department of Veterans Affairs  
DFS/CLIA Certification  
FDA Center for Biologics Evaluation and Research  
FDA Center for Devices and Radiological Health  
FDA Center for Veterinary Medicine  
Health Canada  
Massachusetts Department of Public Health Laboratories  
Ministry of Health and Social Welfare – Tanzania  
National Center of Infectious and Parasitic Diseases (Bulgaria)  
National Health Laboratory Service (South Africa)  
National Institute of Standards and Technology  
National Pathology Accreditation Advisory Council (Australia)  
New York State Department of Health  
Ontario Ministry of Health  
Pennsylvania Dept. of Health  
Saskatchewan Health-Provincial Laboratory  
Scientific Institute of Public Health  
University of Iowa, Hygienic Lab

## Industry Members

3M Medical Division  
AB Biodisk  
Abbott  
Abbott Diabetes Care  
Abbott Molecular Inc.  
Abbott Point of Care Inc.  
Access Genetics  
Acupath  
AdvaMed  
Advancis Pharmaceutical Corporation  
Advantage Bio Consultants, Inc.  
Affymetrix, Inc. (Santa Clara, CA)  
Affymetrix, Inc. (W. Sacramento, CA)  
Agilent Technologies/Molecular Diagnostics  
Ammirati Regulatory Consulting  
Anapharm, Inc.  
Anna Longwell, PC  
Aptium Oncology  
Arpida Ltd.  
A/S Rosco  
Associate Regional & University Pathologists  
Astellas Pharma  
AstraZeneca Pharmaceuticals  
Aviir, Inc.  
Axis-Shield PoC AS  
Bayer Corporation – West Haven, CT  
Bayer HealthCare, LLC, Diagnostics Div. – Elkhart, IN  
BD  
BD Biosciences – San Jose, CA  
BD Diagnostic Systems  
BD Vacutainer Systems  
Beckman Coulter, Inc.  
Beth Goldstein Consultant (PA)  
Bioanalyse, Ltd.  
Bio-Development S.r.l.  
Biomedica Laboratories SDN BHD  
bioMérieux, Inc. (MO)  
bioMérieux, Inc. (NC)  
Bio-Rad Laboratories, Inc. – France  
Bio-Rad Laboratories, Inc. – Irvine, CA  
Bio-Rad Laboratories, Inc. – Plano, TX  
Blaine Healthcare Associates, Inc.  
Braun Biosystems, Inc.  
Canon U.S. Life Sciences, Inc.  
Cempra Pharmaceuticals, Inc.  
Center for Measurement Standards/ITRI  
Centers for Disease Control and Prevention  
Central States Research Centre, Inc.  
Cepheid  
Chen & Chen, LLC (IQUM)  
The Clinical Microbiology Institute  
Comprehensive Cytometric Consulting Control Lab  
Copan Diagnostics Inc.  
Cosmetic Ingredient Review  
Cubist Pharmaceuticals  
Cumbre Inc.  
Dade Behring Marburg GmbH – A  
Siemens Company  
Dahl-Chase Pathology Associates PA  
David G. Rhoads Associates, Inc.  
Diagnostic Products Corporation  
Diagnostica Stago

Dicro, Inc.  
DX Tech  
Eiken Chemical Company, Ltd.  
Elanco Animal Health  
Emisphere Technologies, Inc.  
Eurofins Medinet  
Fio  
Focus Diagnostics  
Future Diagnostics B.V.  
Genomic Health, Inc.  
Gen-Probe  
Genzyme Diagnostics  
GlaxoSmithKline  
GlucoTec, Inc.  
GR Micro LTD  
Greiner Bio-One Inc.  
Habig Regulatory Consulting  
HistoGenex N.V.  
Icon Laboratories, Inc.  
Immunicon Corporation  
Indiana State Department of Health  
Instrumentation Laboratory  
Japan Assn. of Clinical Reagents Industries  
Joanneum Research  
Forschungsgesellschaft mbH  
Johnson & Johnson Pharmaceutical Research and Development, L.L.C.  
Kaiser Permanente  
K.C.J. Enterprises  
Krouwer Consulting  
Laboratory Specialists, Inc.  
LifeScan, Inc. (a Johnson & Johnson Company)  
LipoScience  
Maine Standards Company, LLC  
Medical Device Consultants, Inc.  
Merck & Company, Inc.  
Micromyx, LLC  
MicroPhage  
Monogen, Inc.  
MultiPhase Solutions, Inc.  
Nanogen  
Nanogen, Point-of-Care Diagnostics Div.  
Nanosphere, Inc.  
Nihon Koden Corporation  
Nissui Pharmaceutical Co., Ltd.  
NJK & Associates, Inc.  
NorDx – Scarborough Campus  
NovaBiotics (Aberdeen, UK)  
Novartis Institutes for Biomedical Research  
Nucryst Pharmaceuticals  
Olympus America, Inc.  
Opti Scan Bio Medical Assoc.  
Optimer Pharmaceuticals, Inc.  
Orion Genomics, LLC  
Ortho-Clinical Diagnostics, Inc. (Rochester, NY)  
Ortho-McNeil, Inc.  
Oxonica (UK)  
Panaceapharma Pharmaceuticals  
Paratek Pharmaceuticals, Inc.  
PathCare  
Pathwork Diagnostics  
Pfizer Animal Health  
Pfizer Inc  
Phadia AB  
PlaCor, Inc.  
Powers Consulting Services  
PPD  
ProSource Consulting, Inc.  
QSE Consulting  
Qualtek Clinical Laboratories  
Quest Diagnostics, Incorporated  
Radiometer America, Inc.  
RCC CIDA S. A.  
Replidyne  
Rib-X Pharmaceuticals  
Roche Diagnostics GmbH  
Roche Diagnostics, Inc.  
Roche Diagnostics Ltd  
Roche Diagnostics Shanghai Ltd.  
Roche Molecular Systems  
SAIC Frederick Inc. NCI-Frederick  
Cancer Research & Development Center  
Sanofi Pasteur  
Sarstedt, Inc.  
Schering Corporation  
Sequenom, Inc.  
Siemens Healthcare Diagnostics  
Siemens Medical Solutions Diagnostics (CA)  
Siemens Medical Solutions Diagnostics (DE)  
Siemens Medical Solutions Diagnostics (NY)

Specialty Ranbaxy Ltd  
Sphere Medical Holding Limited  
State of Alabama  
Stirling Medical Innovations  
Streck Laboratories, Inc.  
Sysmex America, Inc. (Mundelein, IL)  
Sysmex Corporation (Japan)  
Targanta Therapeutics, Inc.  
Tethys Bioscience, Inc.  
TheraDoc  
Therapeutic Monitoring Services, LLC  
Theravance Inc.  
Third Wave Technologies, Inc.  
Thrombodyne, Inc.  
ThromboVision, Inc.  
Transasia Bio-Medicals Limited  
Trek Diagnostic Systems  
Upside Endeavors, LLC  
Vital Diagnostics S.r.l.  
Watn-Biolife Diagnostics and Medicals  
Wellstat Diagnostics, LLC  
Wyeth Research  
XDX, Inc.  
YD Consultant  
ZIUR Ltd.

## Trade Associations

AdvaMed  
Japan Association of Clinical Reagents Industries (Tokyo, Japan)

## Associate Active Members

3<sup>rd</sup> Medical Group (AK)  
5th Medical Group/SGSL (ND)  
22 MDSS (KS)  
36th Medical Group/SGSL (Guam)  
48th Medical Group/MDSS (APO, AE)  
55th Medical Group/SGSAL (NE)  
59th MDW/859<sup>th</sup> MDTs/MTL Wilford Hall Medical Center (TX)  
Academisch Ziekenhuis-VUB  
Acadiana Medical Labs, Ltd  
ACL Laboratories (IL)  
Adams County Hospital (OH)  
Air Force Institute for Operational Health (TX)  
Akron's Children's Hospital (OH)  
Alameda County Medical Center  
Albany Medical Center Hospital (NY)  
Albemarle Hospital (NC)  
Alfred I. du Pont Hospital for Children  
All Children's Hospital (FL)  
Allegheny General Hospital (PA)  
Alpena General Hospital (MI)  
Alta Bates Summit Medical Center (CA)  
American University of Beirut Medical Center (NJ)  
Anne Arundel Medical Center (MD)  
Antelope Valley Hospital District (CA)  
Arkansas Children's Hospital (AR)  
Arkansas Dept of Health  
Public Health Laboratory (AR)  
Arkansas Methodist Medical Center (AR)  
Asan Medical Center (Seoul)  
Asante Health System (OR)  
Asiri Group of Hospitals Ltd.  
Asociacion Espanola Primera de Socorros Mutuos (Uruguay)  
Aspirus Wausau Hospital (WI)  
Atlantic City Medical Center (NJ)  
Atlantic Health Sciences Corp.  
Auburn Regional Medical Center (WA)  
Augusta Medical Center (VA)  
Aultman Hospital (OH)  
Avera McKennan (SD)  
Az Sint-Jan  
Azienda Ospedale Di Lecco (Italy)  
Baffin Regional Hospital (Canada)  
Baptist Hospital for Women (TN)  
Baptist Hospital of Miami (FL)  
Basset Army Community Hospital (AK)  
Baton Rouge General (LA)  
Baxter Regional Medical Center (AR)  
Bay Regional Medical Center (MI)  
BayCare Health System (FL)  
Baylor Health Care System (TX)  
Bayou Pathology, APMC (LA)  
Baystate Medical Center (MA)  
B.B.A.G. Ve U. AS., Duzen Laboratories (Turkey)  
Beebe Medical Center (DE)  
Belfast HSS Trust

Beloit Memorial Hospital (WI)	Deaconess Medical Center (WA)	Holy Spirit Hospital (PA)	Maricopa Integrated Health System (AZ)
Ben Taub General Hospital (TX)	Dean Medical Center (WI)	Hopital Cite de La Sante de Laval (Canada)	Marquette General Hospital (MI)
The Bermuda Hospitals Board	DeWitt Healthcare Network (USA)	Hopital du Haut-Richelieu (Canada)	Marshfield Clinic (WI)
Bonnyville Health Center (Canada)	Meddacc (VA)	Hôpital Maisonneuve - Rosemont (Montreal, Canada)	Martha Jefferson Hospital (VA)
Boston Medical Center (MA)	DHHS NC State Lab of Public Health (NC)	Hôpital Sacré-Coeur de Montreal (Quebec, Canada)	Martin Luther King, Jr. Harbor Hospital (CA)
Boulder Community Hospital (CO)	Diagnostic Laboratory Services, Inc. (HI)	Hopital Santa Cabrini Ospedale (Canada)	Martin Memorial Health Systems (FL)
Brantford General Hospital (Canada)	Diagnostic Services of Manitoba (Canada)	Hospital Albert Einstein (Brazil)	Mary Imogene Bassett Hospital (NY)
Bridgeport Hospital (CT)	Diagnósticos da América S/A (Sao Paulo)	Hospital das Clinicas-FMUSP (Brazil)	Marymount Medical Center (KY)
Bronson Methodist Hospital (MI)	DIANON Systems/Lab Corp. (OK)	Hospital de Dirino Espirito Santa (Portugal)	Massachusetts General Hospital (MA)
Broward General Medical Center (FL)	Diaz Gill-Medicina Laboratorial S.A. (MD)	Hôtel Dieu Grace Hospital Library (Windsor, ON, Canada)	Massachusetts General Hospital Division of Laboratory Medicine (MA)
Calgary Laboratory Services (Calgary, AB, Canada)	Dimensions Healthcare System (MD)	Hunter Area Pathology Service (Australia)	Maxwell Air Force Base (AL)
California Pacific Medical Center (CA)	Dr. Erfan & Bagedo General Hospital (Saudi Arabia)	Imelda Hospital (Belgium)	Mayo Clinic (MN)
Cambridge Health Alliance (MA)	Dr. Everette Chalmers Regional Hospital (NB)	Indiana University – Chlamydia Laboratory (IN)	Mayo Clinic Scottsdale (AZ)
Camden Clark Memorial Hospital (WV)	DRAKE Center (OH)	Inova Fairfax Hospital (VA)	MDS Metro Laboratory Services (BC, Canada)
Canadian Science Center for Human and Animal Health (Canada)	Driscoll Children's Hospital (TX)	Institut fur Stand. und Dok. im Med. Lab. (Germany)	Meadows Regional Medical Center (GA)
Cape Breton Healthcare Complex (Canada)	DSI of Bucks County (PA)	Institut National de Santé Publique du Quebec Centre de Doc. – INSPQ (Canada)	Mease Countryside Hospital (FL)
Cape Cod Hospital (MA)	DUHS Clinical Laboratories (NC)	Institute Health Laboratories (PR)	Medecin Microbiologiste (Canada)
Cape Fear Valley Medical Center Laboratory (NC)	Dundee County Hospital (NE)	Institute of Clinical Pathology and Medical Research (Australia)	Medical Center Hospital (TX)
Capital Health/QE II Health Sciences Centre (Nova Scotia)	Durham VA Medical Center (NC)	Institute of Laboratory Medicine Landspítali Univ. Hospital (Iceland)	Medical Center of Louisiana at NO-Charity (LA)
Capital Health - Regional Laboratory Services (Canada)	DVA Laboratory Services (FL)	Institute of Medical & Veterinary Science (Australia)	Medical Centre Ljubljana (Slovenia)
Capital Health System Mercer Campus (NJ)	Dwight D. Eisenhower Medical Center (KS)	Integrated Regional Laboratories South Florida (FL)	Medical College of Virginia Hospital (VA)
Carilion Labs Charlotte	E. A. Conway Medical Center (LA)	International Health Management Associates, Inc. (IL)	Medical Specialists (IN)
Carpenter S.A. de C.V. (Mexico)	East Central Health (Canada)	Ireland Army Community Hospital (KY)	Medical Univ. of South Carolina (SC)
Catholic Health Initiatives (KY)	East Georgia Regional Medical Center (GA)	IWK Health Centre (Canada)	MediCorp - Mary Washington Hospital (VA)
Cavan General Hospital (Ireland)	Eastern Health Pathology (Australia)	Jackson County Memorial Hospital (OK)	Memorial Hermann Healthcare System (TX)
CDC/HIV (APO, AP)	Easton Hospital (PA)	Jackson Health System (FL)	Memorial Hospital at Gulfport (MS)
Cedars-Sinai Medical Center (CA)	Edward Hospital (IL)	Jackson Purchase Medical Center (KY)	Memorial Hospital Laboratory (CO)
Central Baptist Hospital (KY)	Effingham Hospital (GA)	Jacobi Medical Center (NY)	Memorial Medical Center (IL)
Central Kansas Medical Center (KS)	Eliza Coffee Memorial Hospital (AL)	John C. Lincoln Hospital (AZ)	Memorial Medical Center (PA)
Central Texas Veterans Health Care System (TX)	Emory University Hospital (GA)	John Muir Medical Center (CA)	Memorial Regional Hospital (FL)
Centralized Laboratory Services (NY)	Evangelical Community Hospital (PA)	John T. Mather Memorial Hospital (NY)	Mercy Franciscan Mt. Airy (OH)
Centre Hospitalier Anna-Laberge (Canada)	Evans Army Community Hospital (CO)	Johns Hopkins Medical Institutions (MD)	Mercy Hospital (ME)
Centura – Villa Pueblo (CO)	Exeter Hospital (NH)	Johns Hopkins University (MD)	Mercy Medical Center (CO)
Chaleur Regional Hospital (Canada)	Federal Medical Center (MN)	Johnson City Medical Center (TN)	Mercy Medical Center (OR)
Chang Gung Memorial Hospital (Taiwan)	First Health of the Carolinas Moore Regional Hospital (NC)	JPS Health Network (TX)	Methodist Hospital (MN)
Changhua Christian Hospital (Taiwan)	Flaget Memorial Hospital (KY)	Kadlec Medical Center (WA)	Methodist Hospital Pathology (NE)
The Charlotte Hungerford Hospital (CT)	Fletcher Allen Health Care (VT)	Kaiser Permanente (CA)	MetroHealth Medical Center (OH)
Chatham - Kent Health Alliance (Canada)	Fleury S.A. (Brazil)	Kaiser Permanente (MD)	Metropolitan Hospital Center (NY)
Chesapeake General Hospital (VA)	Florida Hospital (FL)	Kaiser Permanente Medical Care (CA)	Metropolitan Medical Laboratory, PLC (IA)
Chester County Hospital (PA)	Florida Hospital Waterman (FL)	Kantonsspital Aarau AG (Switzerland)	The Michener Inst. for Applied Health Sciences (Canada)
Children's Healthcare of Atlanta (GA)	Foot Hospital (MI)	Keller Army Community Hospital (NY)	Middelheim General Hospital
The Children's Hospital (CO)	Fort St. John General Hospital (Canada)	Kenora-Rainy River Reg. Lab. Program (Canada)	Middletown Regional Hospital (OH)
Children's Hospital (OH)	Forum Health Northside Medical Center (OH)	King Fahad National Guard Health Affairs King Abdulaziz Medical City (Saudi Arabia)	Mike O'Callaghan Federal Hospital (NV)
Children's Hospital and Medical Center (WA)	Fox Chase Cancer Center (PA)	King Faisal Specialist Hospital (MD)	Mississippi Baptist Medical Center (MS)
Children's Hospital & Research Center at Oakland (CA)	Frankford Hospital (PA)	King Hussein Cancer Center	Mississippi Public Health Lab (MS)
Children's Hospital Medical Center (OH)	Fraser Health Authority Royal Columbian Hospital Site (Canada)	Kings County Hospital Center (NY)	Monmouth Medical Center (NJ)
Children's Hospital of Philadelphia (PA)	Fresenius Medical Care/Spectra East (NJ)	Kingston General Hospital (Canada)	Montefiore Medical Center (NY)
Children's Hospitals and Clinics (MN)	Fundacio Joan Costa Roma Consorci Sanitari de Terrassa (Spain)	Lab Medico Santa Luzia LTDA (Brazil)	Montreal General Hospital (Quebec)
Children's Medical Center (OH)	Gamma-Dynacare Laboratories (Canada)	Labbette Health (KS)	Morton Plant Hospital (FL)
Children's Medical Center (TX)	Gamma Dynacare Medical Laboratories (Ontario, Canada)	Laboratory Alliance of Central New York (NY)	Mt. Sinai Hospital - New York (NY)
Children's Memorial Hospital (IL)	Garden City Hospital (MI)	LabPlus Auckland Healthcare Services Limited (New Zealand)	Nassau County Medical Center (NY)
The Children's Mercy Hospital (MO)	Garfield Medical Center (CA)	Labway (China)	National Cancer Center (S. Korea)
Childrens Hosp. – Kings Daughters (VA)	Geisinger Medical Center (Danville, PA)	Lafayette General Medical Center (LA)	National Cancer Institute (MD)
Childrens Hospital Los Angeles (CA)	Genesis Healthcare System (OH)	Lakeland Regional Laboratories (MI)	National Healthcare Group (Singapore)
Childrens Hospital of Wisconsin (WI)	George Washington University Hospital (DC)	Lakeland Regional Medical Center (FL)	National Institutes of Health, Clinical Center (MD)
Chilton Memorial Hospital (NJ)	Ghent University Hospital (Belgium)	Lancaster General Hospital (PA)	National Naval Medical Center (MD)
Christiana Care Health Services (DE)	Good Samaritan Hospital (OH)	Landstuhl Regional Medical Center Langley Air Force Base (VA)	National University Hospital Department of Laboratory Medicine (Singapore)
Christus St. John Hospital (TX)	Good Shepherd Medical Center (TX)	LeBonheur Children's Medical Center (TN)	Naval Hospital Great Lakes (IL)
CHU Sainte – Justine (Quebec)	Grana S.A. (Mexico)	Lethbridge Regional Hospital (Canada)	Naval Hospital Oak Harbor (WA)
City of Hope National Medical Center (CA)	Grand Strand Reg. Medical Center (SC)	Lewis-Gale Medical Center (VA)	Naval Medical Center Portsmouth (VA)
Clarian Health – Clarian Pathology Laboratory (IN)	Gundersen Lutheran Medical Center (WI)	L'Hotel-Dieu de Quebec (Quebec, Canada)	NB Department of Health
Cleveland Clinic Health System Eastern Region (OH)	Guthrie Clinic Laboratories (PA)	Licking Memorial Hospital (OH)	The Nebraska Medical Center (NE)
Clinical Labs of Hawaii (HI)	Haga Teaching Hospital (Netherlands)	LifeBridge Health Sinai Hospital (MD)	New England Baptist Hospital (MA)
CLSI Laboratories, Univ. Pittsburgh Med. Ctr. (PA)	Hagerstown Medical Laboratory (MD)	LifeLabs (Canada)	New England Fertility Institute (CT)
Colchester East Hants Health Authority (Canada)	Halton Healthcare Services (Canada)	Loma Linda University Medical (CA)	New Lexington Clinic (KY)
Commonwealth of Virginia (DCLS) (VA)	Hamad Medical Corporation (Qatar)	Long Beach Memorial Medical Center (CA)	New York City Department of Health and Mental Hygiene (NY)
Community Hospital (IN)	Hamilton Regional Laboratory Medicine Program (Canada)	Los Angeles County Public Health Lab. (CA)	New York-Presbyterian Hospital (NY)
The Community Hospital (OH)	Hanover General Hospital (PA)	Louisiana Office of Public Health Laboratory (LA)	New York University Medical Center (NY)
Community Hospital of the Monterey Peninsula (CA)	Harford Memorial Hospital (MD)	Louisiana State University Medical Ctr. (LA)	Newark Beth Israel Medical Center (NJ)
Community Medical Center (NJ)	Harris Methodist Fort Worth (TX)	Lourdes Hospital (KY)	Newton Memorial Hospital (NJ)
Community Memorial Hospital (WI)	Health Network Lab (PA)	Maccabi Medical Care and Health Fund Madison Parish Hospital (LA)	North Bay Hospital (FL)
Consultants Laboratory of W1 LLC (WI)	Health Partners Laboratories Bon Secours Richmond (VA)	Mafraq Hospital	North Carolina Baptist Hospital (NC)
Contra Costa Regional Medical Center (CA)	Health Sciences Research Institute (Japan)	Magnolia Regional Health Center (MS)	North Coast Clinical Laboratory, Inc. (OH)
Cook Children's Medical Center (TX)	Health Waikato (New Zealand)	Main Line Clinical Laboratories, Inc. (PA)	North District Hospital (Hong Kong)
Cork University Hospital (Ireland)	Heartland Health (MO)		North Mississippi Medical Center (MS)
Corpus Christi Medical Center (TX)	Heidelberg Army Hospital (APO, AE)		North Shore-Long Island Jewish Health System Laboratories (NY)
Covance CLS (IN)	Helen Hayes Hospital (NY)		Northeast Pathologists, Inc. (MO)
Covance Evansville (IN)	Hema-Quebec (Canada)		Northridge Hospital Medical Center (CA)
The Credit Valley Hospital (Canada)	Hennepin Faculty Association (MN)		Northside Hospital (GA)
Creighton Medical Laboratories (NE)	Henry Ford Hospital (MI)		Northwest Texas Hospital (TX)
Creighton University Medical Center (NE)	Henry M. Jackson Foundation (MD)		Northwestern Memorial Hospital (IL)
Crozer-Chester Medical Center (PA)	Henry Medical Center, Inc. (GA)		Norton Healthcare (KY)
Darwin Library NT Territory Health Services (Australia)	Hi-Desert Medical Center (CA)		Ochsner Clinic Foundation (LA)
David Grant Medical Center (CA)	Hoag Memorial Hospital Presbyterian (CA)		Ohio State University Hospitals (OH)
Daviness Community Hospital (IN)	Holy Cross Hospital (MD)		Onze Lieve Vrouw Ziekenhuis (Belgium)
Deaconess Hospital Laboratory (IN)	Holy Family Medical Center (WI)		Ordre Professionnel des Technologistes Médicaux du Quebec (Quebec)
	Holy Name Hospital (NJ)		Orebro University Hospital



Orlando Regional Healthcare System (FL)	St. Anthony Hospital (OK)	State of Connecticut Department of Public Health (CT)	University of So. Alabama Children's and Women's Hospital (AL)
The Ottawa Hospital (Canada)	St. Anthony Hospital Central Laboratory (CO)	State of Hawaii Department of Health (HI)	University of Texas Health Center (TX)
Our Lady of Lourdes Medical Center (NJ)	St. Anthony's Hospital (FL)	State of Washington-Public Health Labs (WA)	The University of Texas Medical Branch (TX)
Our Lady of Lourdes Reg. Medical Ctr. (LA)	St. Barnabas Medical Center (NJ)	Steele Memorial Hospital (ID)	University of the Ryukyus (Japan)
Our Lady of the Way Hospital (KY)	St. Christopher's Hospital for Children (PA)	Stillwater Medical Center (OK)	University of Virginia Medical Center
Our Lady's Hospital for Sick Children (Ireland)	St. Elizabeth Community Hospital (CA)	Stony Brook University Hospital (NY)	University of Washington
Overlake Hospital Medical Center (WA)	St. Francis Hospital (SC)	Stormont-Vail Regional Medical Center (KS)	UPMC Bedford Memorial (PA)
Palmetto Health Baptist Laboratory (SC)	St. Francis Medical Center (MN)	Sudbury Regional Hospital (Canada)	U.S.A. Meddae (Pathology Division) (MO)
Pathologists Associated (IN)	St. John Hospital and Medical Center (MI)	Suncoast Medical Clinic (FL)	UW Hospital (WI)
Pathology and Cytology Laboratories, Inc. (KY)	St. John's Hospital (IL)	Sunnybrook Health Science Center (ON, Canada)	UZ-KUL Medical Center (Belgium)
Pathology Associates Medical Laboratories (WA)	St. John's Hospital & Health Ctr. (CA)	Surrise Hospital and Medical Center (NV)	VA (Asheville) Medical Center (NC)
Penn State Hershey Medical Center (PA)	St. John's Mercy Medical Center (MO)	Swedish Medical Center (CO)	VA (Bay Pines) Medical Center (FL)
Pennsylvania Hospital (PA)	St. John's Regional Health Center (MO)	Sydney South West Pathology Service (Australia)	VA (Chillicothe) Medical Center (OH)
Penrose St. Francis Health Services (CO)	St. Joseph Medical Center (MD)	T.J. Samson Community Hospital (KY)	VA (Cincinnati) Medical Center (OH)
The Permanente Medical Group (CA)	St. Joseph Mercy - Oakland (MI)	Taipei Veterans General Hospital (Taiwan)	VA (Dallas) Medical Center (TX)
Perry County Memorial Hospital (IN)	St. Joseph Mercy Hospital (MI)	Taiwan Society of Laboratory Medicine	VA (Dayton) Medical Center (OH)
Peterborough Regional Health Centre (Canada)	St. Joseph's Hospital (FL)	Tallaght Hospital	VA (Decatur) Medical Center (GA)
Piedmont Hospital (GA)	St. Joseph's Hospital & Health Center (ND)	Tartu University Clinics (Tartu)	VA (Hines) Medical Center (IL)
Pitt County Memorial Hospital (NC)	St. Joseph's Medical Center (CA)	Temple Univ. Hospital - Parkinson Pav. (PA)	VA (Indianapolis) Medical Center (IN)
Prairie Lakes Hospital (SD)	St. Joseph's Regional Medical Center (NJ)	Texas Children's Hospital (TX)	VA (Iowa City) Medical Center (IA)
Presbyterian Hospital of Dallas (TX)	St. Jude Children's Research Hospital (TN)	Texas Department of State Health Services	VA (Long Beach) Medical Center (CA)
Presbyterian/St. Luke's Medical Center (CO)	St. Louis University Hospital (MO)	Thomason Hospital (TX)	VA (Miami) Medical Center (FL)
Prince County Hospital	St. Luke's Hospital (FL)	Timmins and District Hospital (Canada)	VA New Jersey Health Care System (NJ)
Princess Margaret Hospital (Hong Kong)	St. Luke's Hospital (IA)	The Toledo Hospital (OH)	VA Outpatient Clinic (OH)
Providence Alaska Medical Center (AK)	St. Luke's Hospital (PA)	Touro Infirmary (LA)	VA (Phoenix) Medical Center (AZ)
Providence Health Care (Canada)	St. Margaret Memorial Hospital (PA)	Tri-Cities Laboratory (WA)	VA (San Diego) Medical Center (CA)
Providence Medford Medical Center (OR)	St. Martha's Regional Hospital (Canada)	Trident Medical Center (SC)	VA (Seattle) Medical Center (WA)
Provincial Health Services Authority (Vancouver, BC, Canada)	St. Mary Medical Center (CA)	Trinity Medical Center (AL)	VA (Sheridan) Medical Center (WY)
Provincial Laboratory for Public Health (Edmonton, AB, Canada)	St. Mary's Health Center (MO)	Tripler Army Medical Center (HI)	VA (Tucson) Medical Center (AZ)
Queen Elizabeth Hospital (Canada)	St. Mary's Healthcare (SD)	Tufts New England Medical Center Hospital (MA)	Valley Health (VA)
Queensland Health Pathology Services (Australia)	St. Mary's Medical Center (IN)	Tulane Medical Center Hospital & Clinic (LA)	Vancouver Hospital and Health Sciences Center (BC, Canada)
Quest Diagnostics, Inc	St. Michael's Hospital Diagnostic Laboratories & Pathology (Canada)	Turku University Central Hospital	Vancouver Island Health Authority (Canada)
Quest Diagnostics, Inc (San Juan Capistrano, CA)	St. Tammany Parish Hospital (LA)	UCI Medical Center (CA)	Vanderbilt University Medical Center (TN)
Quest Diagnostics JV (IN, OH, PA)	St. Thomas More Hospital (CO)	UCLA Medical Center	Via Christi Regional Medical Center (KS)
Quest Diagnostics Laboratories (WA)	Sampson Regional Medical Center (NC)	Clinical Laboratories (CA)	Virga Jesseziekenhuis (Belgium)
Quincy Hospital (MA)	Samsung Medical Center (Korea)	UCSF Medical Center (CA)	ViroMed Laboratories (LabCorp) (MN)
Rady Children's Hospital San Diego (CA)	San Francisco General Hospital-University of California San Francisco (CA)	UMass Memorial Medical Center (MA)	Virtua - West Jersey Hospital (NJ)
Redington-Fairview General Hospital (ME)	Sanford USP Medical Center (SD)	UMC of Southern Nevada (NV)	WakeMed (NC)
Regional Health Authority Four (RHA4) (Canada)	SARL Laboratoire Carron (France)	UNC Hospitals (NC)	Walter Reed Army Medical Center (DC)
Regions Hospital (MN)	Saudi Aramco Medical (Saudi Arabia)	United Clinical Laboratory (Taiwan)	Warren Hospital (NJ)
Reid Hospital & Health Care Services (IN)	Scott Air Force Base (IL)	United Christian Hospital (Hong Kong)	Washington Hospital Center (DC)
Renown Regional Medical Center (NV)	Scott & White Memorial Hospital (TX)	United Clinical Laboratories (IA)	Waterbury Hospital (CT)
Research Medical Center (MO)	Seoul National University Hospital (Korea)	Unity HealthCare (IA)	Waterford Regional Hospital (Ireland)
Riverside Regional Medical Center (VA)	Seton Medical Center (CA)	Universita Campus Bio-Medico (Italy)	Wayne Memorial Hospital (NC)
Riyadh Armed Forces Hospital, Sulaymania	Shamokin Area Community Hospital (PA)	Universitair Ziekenhuis Antwerpen (Belgium)	Weirton Medical Center (WV)
Robert Wood Johnson University Hospital (NJ)	Sheik Kalifa Medical City (UAE)	University College Hospital (Ireland)	Wellstar Douglas Hospital Laboratory (GA)
Roxborough Memorial Hospital (PA)	Shore Memorial Hospital (NJ)	University Medical Center at Princeton (NJ)	Wellstar Paulding Hospital (GA)
Royal Victoria Hospital (Canada)	Shriners Hospitals for Children (SC)	University of Alabama-Birmingham Hospital (AL)	Wellstar Windy Hill Hospital Laboratory (GA)
Rush North Shore Medical Center (IL)	Singapore General Hospital (Singapore)	University of Arkansas for Medical Sci. (AR)	West China Second University Hospital, Sichuan University (P.R. China)
SAAD Specialist Hospital (Saudi Arabia)	SJRMC Plymouth Laboratory (IN)	University of Chicago Hospitals (IL)	West Valley Medical Center Laboratory (ID)
Sacred Heart Hospital (FL)	Sky Lakes Medical Center (OR)	University of Colorado Health Sciences Center (CO)	Westchester Medical Center (NY)
Sacred Heart Hospital (WI)	South Bend Medical Foundation (IN)	University of Colorado Hospital	Western Baptist Hospital (KY)
Sahlgrenska Universitetssjukhuset (Sweden)	South County Hospital (RI)	University of Colorado Health Sciences Center (CO)	Western Healthcare Corporation (Canada)
Saint Elizabeth Regional Medical Center (NE)	South Dakota State Health Laboratory (SD)	University of Iowa Hospitals and Clinics (IA)	Wheaton Franciscan & Midwest Clinical Laboratories (WI)
Saint Francis Hospital & Medical Center (CT)	South Miami Hospital (FL)	University of Kentucky Med. Ctr. (KY)	Wheeling Hospital (WV)
Saint Mary's Regional Medical Center (NV)	Southern Health Care Network (Australia)	University of Maryland Medical System (MD)	Whitehorse General Hospital (Canada)
Saints Memorial Medical Center (MA)	Southern Maine Medical Center (ME)	University of Medicine & Dentistry, NJ University Hosp. (NJ)	William Beaumont Army Medical Center (TX)
St. Agnes Healthcare (MD)	Southwest Nova District Health Authority (Canada)	University of Miami (FL)	William Beaumont Hospital (MI)
	Speare Memorial Hospital (NH)	University of Missouri Hospital (MO)	William Osler Health Center (Canada)
	Spectrum Health - Blodgett Campus (MI)	University of MN Medical Center - Fairview	Winchester Hospital (MA)
	Stanford Hospital and Clinics (CA)	University of MS Medical Center (MS)	Winn Army Community Hospital (GA)
			Wisconsin State Laboratory of Hygiene (WI)
			Wishard Health Sciences (IN)
			Womack Army Medical Center (NC)
			Woodlawn Hospital (IN)
			York Hospital (PA)

## OFFICERS

Robert L. Habig, PhD,  
President  
Habig Regulatory Consulting

Gerald A. Hoeltge, MD,  
President-Elect  
Cleveland Clinic

Mary Lou Gantzer, PhD,  
Secretary  
Dade Behring Inc.

W. Gregory Miller, PhD,  
Treasurer  
Virginia Commonwealth University

Thomas L. Hearn, PhD,  
Immediate Past President  
Centers for Disease Control and Prevention

Glen Fine, MS, MBA,  
Executive Vice President

Susan Blonshine, RRT, RPFT, FAARC  
TechEd

Maria Carballo  
Health Canada

Russel K. Enns, PhD  
Cephaid

Lillian J. Gill, DPA  
FDA Center for Devices and Radiological Health

Prof. Naotaka Hamasaki, MD, PhD  
Nagasaki International University

Valerie Ng, PhD, MD  
Alameda County Medical Center/  
Highland General Hospital

## BOARD OF DIRECTORS

Janet K.A. Nicholson, PhD  
Centers for Disease Control and Prevention

Luann Ochs, MS  
BD Diagnostics - TriPath

Timothy J. O'Leary, MD, PhD  
Department of Veterans Affairs

Klaus E. Stinshoff, Dr.rer.nat.  
Digene (Switzerland) Särl (Retired)

Michael Thein, PhD  
Roche Diagnostics GmbH

James A. Thomas  
ASTM International

# Explore the Latest Offerings From CLSI!

As we continue to set the global standard for quality in laboratory testing, we are adding products and programs to bring even more value to our members and customers.



By becoming a CLSI member, your laboratory will join 1,600+ other influential organizations all working together to further CLSI's efforts to improve health care outcomes. You can play an active role in raising global laboratory testing standards—in your laboratory, and around the world.

Find out which membership option is best for you at [www.clsi.org/membership](http://www.clsi.org/membership).



Find what your laboratory needs to succeed! CLSI U provides convenient, cost-effective continuing education and training resources to help you advance your professional development. We have a variety of easy-to-use, online educational resources that make eLearning stress-free and convenient for you and your staff.

See our current educational offerings at [www.clsi.org/education](http://www.clsi.org/education).



When laboratory testing quality is critical, standards are needed and there is no time to waste. eCLIPSE™ Ultimate Access, our cloud-based online portal of the complete library of CLSI standards, makes it easy to quickly find the CLSI resources you need.

Learn more and purchase eCLIPSE at [cli.org/eCLIPSE](http://cli.org/eCLIPSE).

**For more information, visit [www.clsi.org](http://www.clsi.org) today.**



950 West Valley Road, Suite 2500, Wayne, PA 19087 USA

P: +1.610.688.0100 **Toll Free (US):** 877.447.1888 F: +1.610.688.0700

E: [customerservice@clsi.org](mailto:customerservice@clsi.org) [www.clsi.org](http://www.clsi.org)

ISBN 1-56238-665-4